




DESS

> DENTAL SMART SOLUTIONS



CATALOG 2021/1

TABLE OF CONTENTS

COMPATIBLE COMPONENTS

NOBEL ACTIVE® & NOBEL REPLACE® CC	3
NOBEL REPLACE SELECT™	5
NEODENT® GRAND MORSE	7
STRAUMANN® BONE LEVEL	9
STRAUMANN® BLX	11
STRAUMANN® TISSUE LEVEL & synOcta®	13
3i CERTAIN®	15
ZIMMER SCREW-VENT®	17
BIOHORIZONS® INTERNAL HEX	19
ASTRA EV®	21
ASTRA OSSEOSPEED®	24
MIS® SEVEN INTERNAL HEX	27
MIS® C1	29
MIS® V3	30
MEGAGEN AnyRidge®	31
HIOSEN®	32
CAMLOG®	34
CONELOG®	36
ANKYLOS® C/X	37
DENTIUM®	39
DENTSPLY FRIADENT XIVE®	41

EXTERNAL HEX COMPATIBLE COMPONENTS

NOBEL BRANEMARK®	43
3i OSSEOTITE® EXTERNAL HEX	45
BIOHORIZONS® EXTERNAL HEX	47

FULL MOUTH SOLUTIONS

MULTIUNITS AND COMPONENTS ONLY	48
DESSLOC OVERDENTURE SOLUTIONS	53

CUSTOM MILLING SOLUTIONS

PREMILLED ABUTMENT BLANKS PREMILL	56
ABUTMENT BLANK HOLDERS	59

SCREWS & DRIVER SOLUTIONS

TOOLS, DRIVERS, AND HANDLES	60
RECOMMENDED SCREW TORQUE SETTINGS	61
SCREW SELECTION GUIDE	65

CEREC COMPATIBLE SOLUTIONS

CBASE QUICK REFERENCE- CEREC	67
------------------------------	----

PRODUCT FEATURES

ANALOGS FOR DIGITAL & STONE	68
DESSLOC OVERDENTURE SOLUTIONS	69
LATCH DRIVERS - SUPERIOR DESIGN	70
TIBASE - UNIQUE DESIGN FEATURES	71
SELECT GRIP SURFACE OVERVIEW	72
ANGLEBASE TIBASE SOLUTIONS	73
ELLIPTIBASE FOR NARROW PLATFORMS	74
CBASE- CEREC COMPATIBLE TIBASES	75
PREMILL BLANK HOLDERS	76
MULTIUNIT ABUTMENT SUPERIOR DESIGN	77
NON ENGAGING ANGLE BASE BAR SOLUTIONS	78
LIFETIME WARRANTY	79
PURE SWITCH - 100% CROSS COMPATIBILITY	80
CORPORATE HISTORY	81



ACTIVE HEX

compatible with

Nobel Active® / Nobel Replace® Conical Connection

All abutments include final screw, except all plastic castables in packs of 10 and scan bodys for desktop scanning.

- non-engaging
- engaging

Final Screws - Unigrip		3.0	NP	RP	WP
		SILVER	PURPLE	YELLOW	BLUE
		3.0	3.5	4.3/5.0	5.5
	Single	19.068	19.027	19.030	19.030
	Pack of 10		19.027-P10	19.030-P10	19.030-P10
	Single		19.227	19.230	19.230

Titanium Analog - For 3D Printed Models

	Single	14.068	14.041	14.042	14.069
	Pack of 10		14.041-P10	14.042-P10	

Our analogs are suitable for both 3D and traditional models.

Straight Ti base

	GH				
	1,2mm	●		15.041	15.042
	1,5mm	●	15.068		15.069
	3,0mm	●		15.341	15.342
	1,2mm	●		16.041	16.042
	1,5mm	●	16.068		16.069
	3,0mm	●		16.341	16.342

Ti base Castables

	Pack of 5	●	32.001-P5	32.001-P5	32.001-P5	32.001-P5
		●	33.001-P5	33.001-P5	33.001-P5	33.001-P5

C-Base® - For Cerec

	GH				
	0,8mm	●		58.041	58.042
	3,0mm	●		58.041/3	58.042/3
	0,8mm	●		59.041	59.042
	3,0mm	●		59.041/3	59.042/3

Includes screw, same as standard screw.

AngleBase®		3.0	NP	RP
		SILVER	PURPLE	YELLOW
		3.0	3.5	4.3/5.0
	●		36.041	36.042
	●		37.041	37.042

Plastic sleeves for AngleBase®

	Pack of 5	0°	33.100-P5
	Pack of 5	10°	33.101-P5
	Pack of 5	20°	33.102-P5

ELLIPTIBase®

	Shaft orientation				
	A	Corner	●	56.068	
	B	Flat	●	57.068	

New

Torx® Screw for AngleBase® & ELLIPTIBase®

	(replacement screw)	19.468	19.441	19.442
---	---------------------	--------	--------	--------

Torx® ball Screwdriver	20mm	24mm	30mm	35mm
	DT20TB06	DT24TB06	DT30TB06	DT35TB06

to use with AngleBase® & ELLIPTIBase®

CoCr Abutment

	●		30.041	30.042
	●		31.041	31.042

Includes screw and wax up sleeve

All Plastic Castable Abutments			3.0	NP	RP	WP
			SILVER 3.0	PURPLE 3.5	YELLOW 4.3/5.0	BLUE 5.5
	Single	●		28.041	28.042	
	Pack of 10			11.041-P10	11.042-P10	
	Single	●		29.041	29.042	
	Pack of 10			12.041-P10	12.042-P10	
				Pack of 10 - screw is sold seperately		

Impression Transfer

			17.044	17.045
	Impression transfer packaged with open and closed tray screw and transfer.			

Titanium Temporary Abutments

		●	24.041	24.042
		●	25.041	25.042

Abutments

	Straight		13.041	13.042
	Angled	15°	22.041	22.042
		25°		23.042

Healing Abutments

	3,0mm ^A	20.039/3	20.040/3
	3,0mm ^B	20.047/3	20.048/3
	4,0mm ^A	20.039/4	20.040/4
	4,0mm ^B	20.047/4	20.048/4
	5,0mm ^A	20.039/5	20.040/5
	5,0mm ^B	20.047/5	20.048/5

^Afor single crown restoration

^Bfor bridges and structures restoration

DESSLoc®

	1,0mm	88.039/1	88.040/1
	2,0mm	88.039/2	88.040/2
	3,0mm	88.039/3	88.040/3
	4,0mm	88.039/4	88.040/4

To see the DESSLoc® range please refer to page 53.

Scan Bodys			3.0	NP	RP	WP
			SILVER 3.0	PURPLE 3.5	YELLOW 4.3/5.0	BLUE 5.5
	for desktop scanning * use indicated screw (not included)		50.068 * 19.068	50.041 * 19.027	50.042 * 19.030	50.069 * 19.030
	for intra-oral scanning		52.068	52.041	52.042	52.069
	on Ti base	●	53.100-P10	53.100-P10	53.100-P10	53.100-P10
	Pack of 10	●	53.200-P10	53.200-P10	53.200-P10	53.200-P10

Pre-milled blank Abutments

	Ti	ø10mm	61.068	61.041	61.042	61.069
		ø14mm		62.041	62.042	62.069
	CoCr	ø14mm	71.068	71.041	71.042	71.069

Multi-unit abutments

				MULTI-UNIT PLATFORM 4.8	
	Straight	1,5mm	●	40.007/1	40.008/1
		2,5mm	●	40.007/2	40.008/2
		3,5mm	●	40.007/3	40.008/3
		4,5mm	●		40.008/4
	Angled	17°	●	41.707/2	41.708/2
		3,5mm	●	41.707/3	41.708/3
		30°	●	43.007/3	43.008/3
		4,5mm	●	43.007/4	43.008/4

Angled Multi-units include DLC screw.

To see the components on top of the DESS Multiunits please refer to page 48.



TRI-LOBE

compatible with
Nobel Replace Select™

All abutments include final screw, except castables in packs of 10 and scan bodys for desktop scanning.

Final Screws - Unigrip		NP	RP	WP	6.0
		PURPLE	YELLOW	BLUE	GREEN
	Single	19.004	19.005	19.005	19.005
	Pack of 10	19.004-P10	19.005-P10	19.005-P10	19.005-P10
	DLC coated	19.204	19.205	19.205	19.205

Titanium Analog - For 3D Printed Models

	Single	14.004	14.005	14.006	14.097
	Pack of 10	14.004-P10	14.005-P10		
Our analogs are suitable for both 3D and traditional models.					

Straight Ti base

	GH					
	0,3mm	●	15.004	15.005	15.006	15.097
	1,5mm	●	15.204	15.205	15.206	
	3,0mm	●	15.304	15.305	15.306	
	0,3mm	◀	16.004	16.005	16.006	16.097
	1,5mm	◀	16.204	16.205	16.206	
	3,0mm	◀	16.304	16.305	16.306	

Ti base Castables

	●	32.001-P5	32.001-P5	32.001-P5	32.001-P5
	◀	33.001-P5	33.001-P5	33.001-P5	33.001-P5

C-Base® - For Cerec

	0,4mm	●	58.004		58.006	58.097
	0,5mm	●		58.005		
	1,5mm	●	58.004/2	58.005/2	58.006/2	
	3,0mm	●	58.004/3	58.005/3	58.006/3	
	0,4mm	◀	59.004		59.006	59.097
	0,5mm	◀		59.005		
	1,5mm	◀	59.004/2	59.005/2	59.006/2	
	3,0mm	◀	59.004/3	59.005/3	59.006/3	

Includes screw, same as standard screw.

- non-engaging
- ◀ engaging

AngleBase®

		NP	RP	WP	6.0
		PURPLE	YELLOW	BLUE	GREEN
	●	36.004	36.005	36.006	
	◀	37.004	37.005	37.006	

Plastic sleeves for AngleBase®

Pack of 5	0°	33.100-P5	33.100-P5	33.100-P5
Pack of 5	10°	33.101-P5	33.101-P5	33.101-P5
Pack of 5	20°	33.102-P5	33.102-P5	33.102-P5

the plastic sleeve is non engaging

Torx® Screw for AngleBase®

	(replacement screw)	19.404	19.405	19.405				
				Torx® ball Screwdriver	20mm	24mm	30mm	35mm
					DT20TB06	DT24TB06	DT30TB06	DT35TB06

to use with AngleBase® & ELLIPTIBase®

CoCr Abutment

	●	30.004	30.005	30.006
	◀	31.004	31.005	31.006
Includes screw and wax up sleeve				

New

		NP	RP	WP	6.0
All Plastic Castable Abutments		PURPLE	YELLOW	BLUE	GREEN
	Single	● 28.004	28.005	28.006	
	Pack of 10	11.004-P10	11.005-P10		
	Single	◀ 29.004	29.005	29.006	
	Pack of 10	12.004-P10	12.005-P10		
		Pack of 10 - screw is sold seperately			
Impression Transfer		17.004	17.005	17.006	

 Impression transfer packaged with open and closed tray screw and transfer.

Titanium Temporary Abutments

	● 24.004	24.005	24.006
	◀ 25.004	25.005	25.006

Abutments

	Straight		13.004	13.005	13.006
	Angled	15°	22.004	22.005	22.006
		25°		23.005	

Healing Abutments

	3,0mm	20.004/3	20.005/3	20.006/3
	4,0mm	20.004/4	20.005/4	20.006/4
	5,0mm	20.004/5	20.005/5	20.006/5

DESSLoc®

	1,0mm	88.004/1	88.005/1
	2,0mm	88.004/2	88.005/2
	3,0mm	88.004/3	88.005/3
	4,0mm	88.004/4	88.005/4

To see the DESSLoc® range please refer to page 53.

Scan Bodys

		NP	RP	WP	6.0
		PURPLE	YELLOW	BLUE	GREEN
	for desktop scanning * use indicated screw (not incuded)	50.004 * 19.004	50.005 * 19.005	50.006 * 19.005	50.097 * 19.005
	for intra-oral scanning	52.004	52.005	52.006	52.097
	on Ti base	● 53.100-P10	53.100-P10	53.100-P10	53.100-P10
	Pack of 10	◀ 53.200-P10	53.200-P10	53.200-P10	53.200-P10

Pre-milled blank Abutments

	Ti	ø10mm	61.004	61.005	61.006	61.097
		ø14mm		62.005	62.006	62.097
	CoCr	ø14mm	71.004	71.005	71.006	71.097

Multi-unit abutments

	Straight	1,0mm	●	40.004/1	40.005/1	40.006/1
		2,0mm	●	40.004/2	40.005/2	40.006/2
		3,0mm	●	40.004/3	40.005/3	40.006/3
		4,0mm	●		40.005/4	
		5,0mm	●		40.005/5	
	Angled 17° 	2,0mm	◀	41.704/2	41.705/2	
		3,0mm	◀	41.704/3	41.705/3	
		4,0mm	◀		41.705/4	
	Angled 30° 	4,0mm	◀		43.005/4	
		5,0mm	◀		43.005/5	
		5,0mm	●		43.005/5R	

To see the components on top of the DESS Multiunits please refer to page 48.

All Multi-units include DLC SCREW.



NEO GM



compatible with
Neodent Grand Morse®

All abutments include final screw

Final Screws - .048/1.20 Hex

GM



Single

19.051

Titanium Analog - For 3D Printed Models



14.035

Our analogs are suitable for both 3D and traditional models.

Straight Ti base

GH



0,7mm	●	15.035
1,6mm	●	15.235
3,6mm	●	15.335
0,7mm	●	16.035
1,6mm	●	16.235
3,6mm	●	16.335

Ti base Castables



Pack of 5



32.001-P5
33.001-P5

AngleBase®

GH



2,5mm	●	36.035/3
2,5mm	●	37.035/3

Plastic sleeves for AngleBase®



Pack of 5	0°	33.100-P5
Pack of 5	10°	33.101-P5
Pack of 5	20°	33.102-P5

the plastic sleeve is non engaging

Torx® Screw for AngleBase®



(replacement screw)

19.444

Torx® ball Screwdriver



20mm	24mm	30mm	35mm
DT20TB06	DT24TB06	DT30TB06	DT35TB06

to use with AngleBase® & ELLIPTIBase®

Scan Body



for intra-oral scanning

52.035



CONICAL BL

compatible with
Straumann® Bone level

All abutments include final screw, except all plastic castables in packs of 10 and scan bodys for desktop scanning.

- non-engaging
- engaging

		NC	RC	NC/RC
		3.3	4.1-4.8	NC/RC
	FOR SRA	Single	19.032	19.033
		Pack of 10	19.032-P10	19.037-P10
		Single		19.098
		Pack of 10		19.098-P10

Titanium Analog - For 3D Printed Models		NC	RC
		14.043	14.044
	FOR SRA	Single	
		Pack of 10	14.043-P10
			14.098
Our analogs are suitable for both 3D and traditional models.			

Straight Ti base		GH			
	1,2mm	●	15.043	15.044	
	3,0mm	●	15.343	15.344	
	1,2mm	■	16.043	16.044	
	3,0mm	■	16.343	16.344	
	ON SRA	●			15.098
		■			16.098
	Ti base Castables	●			
		■			
	Pack of 5	●	32.001-P5	32.001-P5	
		■	33.001-P5	33.001-P5	
ON SRA Ti base		●			32.098-P5
		■			33.098-P5

C-Base® - For Cerec				
	0,7mm	●	58.043	58.044
	3,0mm	●	58.043/3	58.044/3
	0,7mm	■	59.043	59.044
	3,0mm	■	59.043/3	59.044/3
Includes screw, same as standard screw.				

AngleBase®		NC	RC
		3.3	4.1-4.8
	●	36.043	36.044
	■	37.043	37.044
Plastic sleeves for AngleBase®			
	Pack of 5	0°	33.100-P5
	Pack of 5	10°	33.101-P5
	Pack of 5	20°	33.102-P5
the plastic sleeve is non engaging			

ELLIPTIBase®		Shaft orientation	
	A Corner	✓	56.043
	B Flat	✗	57.043

New


Torx® Screw for AngleBase® & ELLIPTIBase®			
	(replacement screw)	19.443	19.443
Torx® ball Screwdriver		20mm	24mm
		DT20TB06	DT24TB06
		DT30TB06	DT35TB06
to use with AngleBase® & ELLIPTIBase®			

CoCr Abutment			
	Includes screw and wax up sleeve	●	30.043
		■	31.043
ON ABUTMENT SRA			
			30.098
	Includes screw and wax up sleeve		31.098

New

		NC	RC
All Plastic Castable Abutments		3.3	4.1-4.8
	Single	●	28.043 28.044
	Pack of 10		11.043-P10 11.044-P10
	Single	●	29.043 29.044
	Pack of 10		12.043-P10 12.044-P10

Impression Transfer

		17.046	17.047
Impression transfer packaged with open and closed tray screw and transfer.			

Titanium Temporary Abutments

		●	24.043 24.044
		●	25.043 25.044
		●	24.098
		●	25.098

Abutments

	Straight		13.043	13.044
	Angled	15°	22.043	22.044

Healing Abutments

	3,0mm	20.041/3	20.042/3
	4,0mm	20.041/4	20.042/4
	5,0mm	20.041/5	20.042/5

DESSLoc®

	1,0mm		88.042/1
	2,0mm	88.041/2	88.042/2
	3,0mm	88.041/3	88.042/3
	4,0mm	88.041/4	88.042/4

To see the DESSLoc® range please refer to page 53.

Scan Bodys

		NC	RC
		3.3	4.1-4.8
	for desktop scanning * use indicated screw (not incuded)	50.043 * 19.032	50.044 * 19.033
	for intra-oral scanning	52.043	52.044
	ON ABUTMENT SRA for desktop scanning * use indicated screw (not incuded)		50.098 * 19.098
	for intra-oral scanning		52.098
	on Ti base		
	Pack of 10	●	53.100-P10 53.100-P10
		●	53.200-P10 53.200-P10

Pre-milled blank Abutments

	Ti	ø10mm	61.043	61.044
		ø14mm	62.043	62.044
	CoCr	ø14mm	71.043	71.044

To see all the Pre-milled blank available please refer to page 138.

Multi-unit abutments

					
Straight	1,5mm		40.043/1	40.044/1	
	2,5mm		40.043/2	40.044/2	
	3,5mm		40.043/3	40.044/3	
	4,5mm			40.044/4	
					
Angled 17°	2,5mm		41.743/2	41.744/2	
	3,5mm		41.743/3	41.744/3	
					
30°	3,5mm		43.043/3	43.044/3	
	4,5mm		43.043/4	43.044/4	

To see the components on top of the DESS Multiunits please refer to page 48.



CONICAL BLX

compatible with
Straumann® BLX



RB WB RB/WB

Final Screws - Torx®

PURPLE BROWN



19.054	19.054	19.054
--------	--------	--------

Titanium Analog - For 3D Printed Models



14.062	14.063
--------	--------

Our analogs are suitable for both 3D and traditional models.

Straight Ti base



0,8mm	●		15.063	
1,5mm	●	15.062/2	15.063/2	15.062/2
2,5mm	●	15.062/3		15.062/3
3,5mm	●	15.062/4		15.062/4
0,8mm	★		16.063	
1,5mm	★	16.062/2	16.063/2	16.062/2
2,5mm	★	16.062/3		16.062/3
3,5mm	★	16.062/4		16.062/4

Ti base Castables



Pack of 5	●	32.001-P5
	★	33.001-P5

RB WB RB/WB

PURPLE BROWN

AngleBase®



36.062	36.063	36.062
37.062	37.063	37.062

Plastic sleeves for AngleBase®



Pack of 5	0°	33.100-P5
Pack of 5	10°	33.101-P5
Pack of 5	20°	33.102-P5

the plastic sleeve is non engaging

Torx® Screw for AngleBase®



19.441	19.441	19.441
--------	--------	--------

Torx® ball Screwdriver	20mm	24mm	30mm	35mm
	DT20TB06	DT24TB06	DT30TB06	DT35TB06

to use with AngleBase® & ELLIPTIBase®

Healing Abutments



0,8mm ^{A/B}		20.463/1	
1,5mm ^A	20.461/2	20.461/2	20.461/2
1,5mm ^B	20.462/2	20.462/2	20.462/2
1,5mm ^{A/B}		20.463/2	
2,5mm ^A	20.461/3	20.461/3	20.461/3
2,5mm ^B	20.462/3	20.462/3	20.462/3
3,5mm ^A	20.461/4	20.461/4	20.461/4
3,5mm ^B	20.462/4	20.462/4	20.462/4

^A for single crown restoration

^B for bridges and structures restoration




OCTAGON

compatible with
Straumann® Tissue level & synOcta®

All abutments include final screw, except all plastic castables in packs of 10 and scan bodys for desktop scanning.

- non-engaging
- engaging

		NNC	RN	WN
Final Screws - Torx®		3.3	4.1/4.8	4.8/6.5
	Single	19.070	19.026	19.026
	Pack of 10		19.026-P10	19.026-P10
ON ABUTMENT				
	Single		19.009	19.009
	Pack of 10		19.009-P10	19.009-P10
Titanium Analog - For 3D Printed Models				
	Single	14.070	14.009	14.045
	Pack of 10		14.009-P10	
ON ABUTMENT				
	Single		14.010	14.046
	Pack of 10		14.010-P10	
Straight Ti base		GH		
	0,3mm	●		15.045
	0,4mm	●		15.046
	0,5mm	●	15.070	
	0,3mm	●		16.045
	0,4mm	●		16.009
	0,5mm	●	16.070	
ON ABUTMENT				
	0,3mm	●	15.010	15.048
	0,3mm	●	16.010	16.048
Ti base Castables				
	Pack of 5	●	32.046-P5	32.045-P5
		●	33.009-P5	33.045-P5
ON ABUTMENT				
	Pack of 5	●	32.010-P5	32.045-P5
		●	33.010-P5	33.045-P5

		NNC	RN	WN
C-Base® - For Cerec		3.3	4.1/4.8	4.8/6.5
	1,0mm	●	58.009	
	1,2mm	●		58.045
	1,0mm	●	59.009	
	1,2mm	●		59.045
Includes screw, same as standard screw.				

		●		
AngleBase®				
		●	36.046	36.047
		●	37.009	37.045
Plastic sleeves for AngleBase®				
Pack of 5	0°		33.200-P5	
Pack of 5	10°		33.201-P5	
Pack of 5	20°		33.301-P5	
the plastic sleeve is non engaging				

		Shaft orientation		
ELLIPTIBase®				
	A Corner		56.070	
	B Flat		57.070	

Torx® Screw for AngleBase® & ELLIPTIBase®						
	(replacement screw)		19.446	19.446		
Torx® ball Screwdriver		20mm	24mm	30mm	35mm	
		DT20TB06	DT24TB06	DT30TB06	DT35TB06	
to use with AngleBase® & ELLIPTIBase®						

ON ABUTMENT - connects onto the top of the Straumann 1.5 or DESS Converter Abutment (on page 14). These abutments are similar but not exactly like a straight multiunit. Please see the following pages for parts that will fit on top of these abutments.

		RN	WN
CoCr Abutment		4.1/4.8	4.8/6.5
	ON ABUTMENT	30.009	30.045
		31.009	31.045
	Includes screw and wax up sleeve	30.010	30.048
		31.010	31.048

All Plastic Castable Abutments

	Single	●	28.009	28.045
	Pack of 10		11.009-P10	
	Single	●	29.009	29.045
	Pack of 10		12.009-P10	
ON ABUTMENT				
	Single	●	28.010	28.046
	Pack of 10		11.010-P10	
	Single	●	29.010	29.046
	Pack of 10		12.010-P10	

Impression Transfer

			17.009	17.053
	ON ABUTMENT			
			17.010	17.048
Impression transfer packaged with open and closed tray screw and transfer.				

Titanium Temporary Abutments

	ON ABUTMENT		24.026	24.045
			25.009	25.045
			24.010	24.047
			25.026	25.047

Abutments

	Straight		13.009*	
	Angled	15°	22.009*	
		25°	23.009*	
			*inc.captive screw	

Healing Abutments

	2,0mm	20.009/2	20.045/2
	3,0mm	20.009/3	20.045/3
	4,0mm	20.009/4	20.045/4
	ON ABUTMENT		20.010

Converter Abutment

	ON ABUTMENT		35.001	
			34.001	34.002

DESSLoc®

	1,0mm	88.009/1	
	2,0mm	88.009/2	
	3,0mm	88.009/3	
	4,0mm	88.009/4	
To see the DESSLoc® range please refer to page 53.			

Scan Bodys

	for desktop scanning	50.070	50.009	50.045
	* use indicated screw not included	* 19.070	* 19.026	* 19.026
	ON ABUTMENT			
			50.010	50.048
for intra-oral scanning			* 19.009	* 19.009
ON ABUTMENT		52.070	52.009	52.045
on Ti base			52.010	52.048
Pack of 10			53.100-P10	
			53.200-P10	

Pre-milled blank Abutments

	Ti	ø10mm	61.070	61.009	61.045
		ø14mm		62.009	62.045
	CoCr	ø14mm	71.070	71.009	71.045

ON ABUTMENT - connects onto the top of the Straumann 1.5 or DESS Converter Abutment (on this page). This abutment is similar to a straight multiunit abutment. Please see the following pages for parts that will fit on top of these abutments.



INTERNAL HEX "CLICK"

compatible with
3i Certain®

All abutments include final screw, except all plastic castables in packs of 10 and scan bodys for desktop scanning.

- non-engaging
- engaging

		3.4	4.1	5.0
		PURPLE	BLUE	YELLOW
	Single	19.011	19.011	19.011
	Pack of 10	19.011-P10	19.011-P10	19.011-P10

Titanium Analog - For 3D Printed Models

	Single	14.014	14.015	14.016
	Pack of 10		14.015-P10	
Our analogs are suitable for both 3D and traditional models.				

Straight Ti base

	0,3mm	●	15.014	15.015	15.016
	1,5mm	●	15.214	15.215	15.216
	3,0mm	●	15.314	15.315	15.316
	0,3mm	●	16.014	16.015	16.016
	1,5mm	●	16.214	16.215	16.216
	3,0mm	●	16.314	16.315	16.316

Ti base Castables



●	32.001-P5	32.001-P5	32.001-P5
●	33.001-P5	33.001-P5	33.001-P5

C-Base® - For Cerec



0,4mm	●	58.014	58.015	
0,5mm	●			58.016
1,5mm	●	58.014/2	58.015/2	58.016/2
3,0mm	●	58.014/3	58.015/3	58.016/3
0,4mm	●	59.014	59.015	
0,5mm	●			59.016
1,5mm	●	59.014/2	59.015/2	59.016/2
3,0mm	●	59.014/3	59.015/3	59.016/3

Includes screw, same as standard screw.

AngleBase®



3.4 4.1 5.0
PURPLE BLUE YELLOW

●	36.014	36.015	36.016
●	37.014	37.015	37.016

Plastic sleeves for AngleBase®



Pack of 5	0°	33.100-P5	33.100-P5	33.100-P5
Pack of 5	10°	33.101-P5	33.101-P5	33.101-P5
Pack of 5	20°	33.102-P5	33.102-P5	33.102-P5

the plastic sleeve is non engaging

Torx® Screw for AngleBase®



19.438	19.438	19.438
--------	--------	--------

Torx® ball Screwdriver 20mm 24mm 30mm 35mm



DT20TB06	DT24TB06	DT30TB06	DT35TB06
----------	----------	----------	----------

to use with AngleBase® & ELLIPTIBase®

CoCr Abutment




●	30.014	30.015	30.016
●	31.014	31.015	31.016

Includes screw and wax up sleeve

New

All Plastic Castable Abutments			3.4	4.1	5.0
			PURPLE	BLUE	YELLOW
	Single	●	28.014	28.015	28.016
				11.015-P10	
	Pack of 10	●	29.014	29.015	29.016
				12.015-P10	

Impression Transfer

		17.014	17.015	17.016
	Impression transfer packaged with open and closed tray screw and transfer.			

Titanium Temporary Abutments

		●	24.014	24.015	24.016
		●	25.014	25.015	25.016

Abutments

	Straight		13.014	13.015	13.016
	Angled	15°	22.014	22.015	22.016
		25°		23.015	

Healing Abutments			3.4	4.1	5.0
			PURPLE	BLUE	YELLOW
	3,0mm		20.014/3	20.015/3	20.016/3
	4,0mm		20.014/4	20.015/4	20.016/4
	5,0mm		20.014/5	20.015/5	20.016/5

DESSLoc®

	1,0mm		88.014/1	88.015/1	
	2,0mm		88.014/2	88.015/2	
	3,0mm		88.014/3	88.015/3	
	4,0mm		88.014/4	88.015/4	
To see the DESSLoc® range please refer to page 53.					

Scan Bodys

	for desktop scanning * use indicated screw not included		50.014 * 19.011	50.015 * 19.011	50.016 * 19.011
	for intra-oral scanning		52.014	52.015	52.016
	on Ti base	Pack of 10	●	53.100-P10	53.100-P10
			●	53.200-P10	53.200-P10

Pre-milled blank Abutments

	Ti	ø10mm	61.014	61.015	61.016
		ø14mm		62.015	62.016
	CoCr	ø14mm	71.014	71.015	71.016



INTERNAL HEX USA

compatible with
Zimmer® Screw-Vent

- non-engaging
- engaging

All abutments include final screw, except all plastic castables in packs of 10 and scan bodys for desktop scanning.

Final Screws - .050/1.27Hex

GREEN PURPLE YELLOW

3.5 4.5 5.7



Single	19.042	19.042	19.042
Pack of 10	19.042-P10	19.042-P10	19.042-P10

Titanium Analog - For 3D Printed Models



Single	14.017	14.018	14.019
Pack of 10	14.017-P10	14.018-P10	

Our analogs are suitable for both 3D and traditional models.

Straight Ti base



GH				
0,3mm	●		15.018	15.019
0,8mm	●	15.017		
1,5mm	●	15.217	15.218	15.219
3,0mm	●	15.317	15.318	15.319
0,3mm	●		16.018	16.019
0,8mm	●	16.017		
1,5mm	●	16.217	16.218	16.219
3,0mm	●	16.317	16.318	16.319

Ti base Castables



Pack of 5	●	32.001-P5	32.001-P5	32.001-P5
	●	33.001-P5	33.001-P5	33.001-P5

C-Base® - For Cerec

GH				
0,5mm	●	58.017	58.018	58.019
1,5mm	●	58.017/2	58.018/2	58.019/2
3,0mm	●	58.017/3	58.018/3	58.019/3
0,5mm	●	59.017	59.018	59.019
1,5mm	●	59.017/2	59.018/2	59.019/2
	●		59.018/3	59.019/3

New



AngleBase®



●	36.017	36.018	36.019
●	37.017	37.018	37.019


Plastic sleeves for AngleBase®



Pack of 5	0°	33.100-P5
Pack of 5	10°	33.101-P5
Pack of 5	20°	33.102-P5

the plastic sleeve is non engaging

ELLIPTIBase®



Shaft orientation		
A Corner	●	56.017
B Flat	●	57.017

New

Torx® Screw for AngleBase® & ELLIPTIBase®



(replacement screw)

19.417	19.417	19.417
--------	--------	--------

Torx® ball Screwdriver	20mm	24mm	30mm	35mm
	DT20TB06	DT24TB06	DT30TB06	DT35TB06

to use with AngleBase® & ELLIPTIBase®

CoCr Abutment



●	30.017	30.018	30.019
●	31.017	31.018	31.019

Includes screw and wax up sleeve

			GREEN	PURPLE	YELLOW
All Plastic Castable Abutments			3.5	4.5	5.7
	Single	●	28.017	28.018	28.019
	Pack of 10		11.017-P10	11.018-P10	
	Single	●	29.017	29.018	29.019
	Pack of 10		12.017-P10	12.018-P10	

Impression Transfer

			17.017	17.018	17.019
	Impression transfer packaged with open and closed tray screw and transfer.				

Titanium Temporary Abutments

		●	24.017	24.018	24.019
		●	25.017	25.018	25.019

Abutments

	Straight		13.017	13.018	13.019
	Angled	15°	22.017	22.018	22.019
		25°	23.017	23.018	

Healing Abutments

	3,0mm	20.017/3	20.018/3	20.019/3
	4,0mm	20.017/4	20.018/4	20.019/4
	5,0mm	20.017/5	20.018/5	20.019/5

DESSLoc®

	1,0mm	88.017/1	88.018/1	
	2,5mm	88.017/2	88.018/2	
	3,5mm	88.017/3	88.018/3	
	4,5mm	88.017/4	88.018/4	
	5,5mm	88.017/5		

To see the DESSLoc® range please refer to page 53.

			GREEN	PURPLE	YELLOW
Scan Bodys			3.5	4.5	5.7
	for desktop scanning * use indicated screw not included		50.017 * 19.012	50.018 * 19.012	50.019 * 19.012
	for intra-oral scanning		52.017	52.018	52.019
	on Ti base Pack of 10	●	53.100-P10	53.100-P10	53.100-P10
		●	53.200-P10	53.200-P10	53.200-P10

Pre-milled blank Abutments

	Ti	ø10mm	61.017	61.018	61.019
		ø14mm	62.017	62.018	62.019
	CoCr	ø14mm	71.017	71.018	71.019
			MULTI-UNIT PLATFORM		

Multi-unit abutments

				4.8	
	Straight	1,0mm	●	40.019/1	40.020/1
		2,0mm	●	40.019/2	40.020/2
		3,0mm	●	40.019/3	40.020/3
		4,0mm	●	40.019/4	
		5,0mm	●	40.019/5	
	Angled 17°	2,5mm	●	41.719/2	41.720/2
		3,5mm	●	41.719/3	41.720/3
	30°	3,5mm	●	43.019/3	43.020/3
		4,5mm	●	43.019/4	43.020/4

Angled Multi-units include DLC screw.

To see the components on top of the DESS Multiunits please refer to page 48.



INTERNAL HEX BH

compatible with

BioHorizons® Internal Hex


- non-engaging
- ⬡ engaging

Final Screws - .050/1.27Hex



	SILVER 3.0	YELLOW 3.5	GREEN 4.5	BLUE 5.7
Single	19.043	19.043	19.043	19.043

Titanium Analog - For 3D Printed Models



Single	14.099	14.087	14.088	14.089
--------	--------	--------	--------	--------

Our analogs are suitable for both 3D and traditional models.

Straight Ti base



GH				
0,3mm	●		15.088	15.089
0,8mm	●		15.087	
0,3mm	⬡		16.088	16.089
0,8mm	⬡		16.087	

Ti base Castables



Pack of 5	●			
	⬡	32.001-P5	32.001-P5	32.001-P5
		33.001-P5	33.001-P5	33.001-P5

C-Base® - For Cerec



0,5mm	●	58.087	58.088	58.089
3,0mm	●	58.087/3	58.088/3	58.089/3
0,5mm	⬡	59.087	59.088	59.089
3,0mm	⬡	59.087/3	59.088/3	59.089/3

Includes screw, same as standard screw.

AngleBase®



SILVER 3.0	YELLOW 3.5	GREEN 4.5	BLUE 5.7
---------------	---------------	--------------	-------------

●	36.087	36.088	36.089
⬡	37.087	37.088	37.089

Plastic sleeves for AngleBase®



Pack of 5	0°	33.100-P5	33.100-P5
Pack of 5	10°	33.101-P5	33.101-P5
Pack of 5	20°	33.102-P5	33.102-P5

ELLIPTIBase®



Shaft orientation

A Corner	⬡	56.099
B Flat	⬡	57.099

New

Torx® Screw for AngleBase®



(replacement screw)

Torx® ball Screwdriver



19.417	19.417	19.417	19.417
20mm	24mm	30mm	35mm
DT20TB06	DT24TB06	DT30TB06	DT35TB06

to be used with AngleBase® & ELLIPTIBase®

Healing Abutments



3,0mm	20.017/3	20.018/3	20.019/3
4,0mm	20.017/4	20.018/4	20.019/4
5,0mm	20.017/5	20.018/5	20.019/5

New

DESSLoc®		YELLOW	GREEN
		3.5	4.5
	1,0mm	88.017/1	88.018/1
	2,5mm	88.017/2	88.018/2
	3,5mm	88.017/3	88.018/3
	4,5mm	88.017/4	88.018/4
	5,5mm	88.017/5	

To see the DESSLoc® range please refer to page 53.

Multi-unit abutments		MULTI-UNIT PLATFORM			
		4.8			
	Straight	1,0mm	●	40.019/1	40.020/1
		2,0mm	●	40.019/2	40.020/2
		3,0mm	●	40.019/3	40.020/3
		4,0mm	●	40.019/4	
		5,0mm	●	40.019/5	
	Angled	2,5mm	⬡	41.719/2	41.720/2
		3,5mm	⬡	41.719/3	41.720/3
		3,5mm	⬡	43.019/3	43.720/3
		4,5mm	⬡	43.019/4	43.720/4

To see the components on top of the DESS Multiunits please refer to page 48.

Scan Bodys		SILVER	YELLOW	GREEN	BLUE
		3.0	3.5	4.5	5.7
	for intra-oral scanning	52.099	52.087	52.088	52.089
	on Ti base				
	Pack of 10	●	53.100-P10	53.100-P10	53.100-P10
		⬡	53.200-P10	53.200-P10	53.200-P10

Pre-milled blank Abutments		Ti	ø10mm				
			ø14mm				
				61.099	61.087	61.088	61.089
						62.088	62.089

INTERNAL HEX BH



CONIC EVO

compatible with
Astra® EV

All abutments include final screw, except scan bodys for desktop scanning.

Final Screws - .050/1.27Hex		3.6 PURPLE	4.2 YELLOW	4.8 BLUE
	Single	19.059	19.060	19.061
	ON UNIABUTMENT®	19.058		

Titanium Analog - For 3D Printed Models		3.6 PURPLE	4.2 YELLOW	4.8 BLUE
	Single	14.059	14.060	14.061
	OF UNIABUTMENT®	14.058		
		Our analogs are suitable for both 3D and traditional models.		



1,0mm	●	15.059	15.060	15.061
3,0mm	●	15.359	15.360	15.361
1,0mm	⚙	16.059	16.060	16.061
3,0mm	⚙	16.359	16.360	16.361
ON UNIABUTMENT®				
	●	15.058		



Pack of 5	●	32.001-P5	32.001-P5	32.001-P5
	⚙	33.001-P5	33.001-P5	33.001-P5
ON UNIABUTMENT®				
Pack of 5	●	32.001-P5		

C-Base® - For Cerec



1,0mm		58.059	58.060	58.061
3,0mm		58.059/3	58.060/3	58.061/3
1,0mm		59.059	59.060	59.061
3,0mm		59.059/3	59.060/3	59.061/3

Includes screw, same as standard screw.

- non-engaging
- ⚙ engaging

AngleBase®		3.6 PURPLE	4.2 YELLOW	4.8 BLUE
	●	36.059	36.060	36.061
	⚙	37.059	37.060	37.061
ON UNIABUTMENT®				
		●	36.058	

New

Plastic sleeves for AngleBase®				
Pack of 5	0°	33.100-P5	33.100-P5	33.100-P5
Pack of 5	10°	33.101-P5	33.101-P5	33.101-P5
Pack of 5	20°	33.102-P5	33.102-P5	33.102-P5

Torx® Screw for AngleBase®				
	(replacement screw)	19.441	19.460	19.461
	ON UNIABUTMENT®			
		(replacement screw)	19.331	

New

Torx® ball Screwdriver	20mm	24mm	30mm	35mm
	DT20TB06	DT24TB06	DT30TB06	DT35TB06

to use with AngleBase® & ELLIPTIBase®

CoCr Abutment				
	●	30.059	30.060	30.061
	⬢	31.059	31.060	31.061
	ON UNIABUTMENT®			
	●	30.058		
Includes screw and wax up sleeve				

All Plastic Castable Abutments		3.6 PURPLE	4.2 YELLOW	4.8 BLUE
	Single 	28.059	28.060	28.061
	ON UNIABUTMENT®			
	Single  Pack of 10	28.058		

Impression Transfer

		17.059	17.060	17.061
	ON UNIABUTMENT®			
		17.058		

Impression transfer packaged with open and closed tray screw and transfer.

Titanium Temporary Abutments

		24.059	24.060	24.061
		25.059	25.060	25.061
	ON UNIABUTMENT®			
	 24.058			

Healing Abutments

	3,5mm	ø4 20.059/3	ø5 20.060/3	ø5 20.061/3
	4,5mm	ø4 20.059/4	ø5 20.060/4	ø6,5 20.061/4
	6,5mm	ø4 20.059/6	ø5 20.060/6	ø6,5 20.061/6

DESSLoc®		3.6 PURPLE	4.2 YELLOW	4.8 BLUE
	1,0mm	88.059/1	88.060/1	88.061/1
	2,0mm	88.059/2	88.060/2	88.061/2
	3,0mm	88.059/3	88.060/3	88.061/3
	4,0mm	88.059/4	88.060/4	88.061/4

To see the DESSLoc® range please refer to page 53.

Scan Bodys

	for desktop scanning * use indicated screw (not included)	50.059 * 19.059	50.060 * 19.060	50.061 * 19.061
	ON UNIABUTMENT®			
		50.058 * 19.058		
	for intra-oral scanning	52.059	52.060	52.061
	ON UNIABUTMENT®			
		52.058		
	on Ti base			
	Pack of 10 	53.100-P10	53.100-P10	53.100-P10
		53.200-P10	53.200-P10	53.200-P10

Pre-milled blank Abutments

	Ti	ø10mm	61.059	61.060	61.061
		ø14mm	62.059	62.060	62.061
	CoCr	ø14mm	71.059	71.060	71.061

UniAbutment®

3.6
PURPLE

4.2
YELLOW

4.8
BLUE

	UniAbutment needs a 1,75mm hex screwdriver	1,0mm	●	40.059/1	40.060/1	40.061/1
		2,0mm	●	40.059/2	40.060/2	40.061/2
		3,0mm	●	40.059/3	40.060/3	40.061/3
		5,0mm	●	40.059/5	40.060/5	40.061/5

Multi-unit abutments

	Straight		1,5mm	●	40.027/1	40.036/1	40.037/1
			2,5mm	●	40.027/2	40.036/2	40.037/2
			3,5mm	●	40.027/3	40.036/3	40.037/3
	Angled	17°	3,5mm	⚙	41.727/3	41.736/3	41.737/3
			4,5mm	⚙	41.727/4	41.736/4	41.737/4
		30°	4,5mm	⚙	43.027/4	43.036/4	43.037/4
			5,5mm	⚙	43.027/5	43.036/5	43.037/5

To see the components on top of the DESS Multiunits please refer to page 48.



INTERNAL HEX CONIC

compatible with
Astra Tech Osseospeed™

All abutments include final screw, except all plastic castables in packs of 10 and scan bodys for desktop scanning.

- non-engaging
- engaging

Final Screws - .050/1.27Hex

YELLOW 3.0AQUA 3.5-4.0LILAC 4.5-5.0

	Single	19.015	19.016	19.017
	Pack of 10		19.016-P10	19.017-P10
ON UNIABUTMENT®				
	Single	19.018		
	Pack of 10	19.018-P10		

Titanium Analog - For 3D Printed Models

	OF UNIABUTMENT®	Single	14.023	14.024	14.025
		Pack of 10		14.024-P10	14.025-P10
		20° Single	14.036		
		Pack of 10	14.036-P10		
		45° Single	14.037		
		Pack of 10	14.037-P10		

Straight Ti base

	0,9mm	●		15.025
	1,0mm	●	15.023	
	1,1mm	●		15.024
	3,0mm	●	15.324	15.325
	0,9mm	○		16.025
	1,0mm	○	16.023	
	1,1mm	○		16.024
	3,0mm	○	16.324	16.325
	ON UNIABUTMENT®			
	20°	●	15.036	
	45°	●	15.037	

Ti base Castables

		32.001-P5	32.001-P5	32.001-P5
		33.001-P5	33.001-P5	33.001-P5
	ON UNIABUTMENT®	32.001-P5		

C-Base® - For Cerec

AQUA 3.5-4.0LILAC 4.5-5.0

	0,8mm	●	58.024	
	0,9mm	●		58.025
	3,0mm	●	58.024/3	58.025/3
	0,8mm	○	59.024	
	0,9mm	○		59.025
	3,0mm	○	59.024/3	59.025/3
Includes screw, same as standard screw.				

New

AngleBase®

	●	36.024	36.025
	○	37.024	37.025

Plastic sleeves for AngleBase®

	Pack of 5	0°	33.100-P5	33.100-P5	33.100-P5
	Pack of 5	10°	33.101-P5	33.101-P5	33.101-P5
	Pack of 5	20°	33.102-P5	33.102-P5	33.102-P5

Torx® Screw for AngleBase®


	(replacement screw)				19.441	19.425
Torx® ball Screwdriver		20mm	24mm	30mm	35mm	
		DT20TB06	DT24TB06	DT30TB06	DT35TB06	

CoCr Abutment to use with AngleBase® & ELLIPTIBase®

	●	30.024	30.025
	○	31.024	31.025
ON UNIABUTMENT®			
	20°	●	30.036
	45°	●	30.037

Includes screw and wax up sleeve

All Plastic Castable Abutments		AQUA 3.5-4.0		LILAC 4.5-5.0
	Single	●	28.024	28.025
	Pack of 10		11.024-P10	11.025-P10
	Single	○	29.024	29.025
	Pack of 10		12.024-P10	12.025-P10
	ON UNIABUTMENT®			
	20° Single	●	28.036	
45°	Pack of 10		11.036-P10	
	Single	●	28.037	
	Pack of 10		11.037-P10	

Impression Transfer					
	ON UNIABUTMENT®		20°	17.036	17.024
			45°	17.037	17.025
	Impression transfer packaged with open and closed tray screw and transfer.				

Titanium Temporary Abutments					
			●	24.024	24.025
			○	25.024	25.025
	ON UNIABUTMENT®				
	20°	●	24.028		
45°	●	24.029			

Abutments					
	Straight		13.024	13.025	
	Angled	15°	22.024	22.025	
		25°	23.024		

Healing Abutments		YELLOW 3.0		AQUA 3.5-4.0	LILAC 4.5-5.0
	Anatomic Healing abutment	3,0mm		20.024/3	20.025/3
		4,0mm		20.024/4	20.025/4
		5,0mm		20.024/5	20.025/5
		4,0mm		Ø4,5 20.424/4	
		4,0mm		Ø5,5 20.524/4	20.525/4
		4,0mm			Ø6,5 20.625/4
	DESSLoc®				
	1mm	88.023/1	88.024/1	88.025/1	
	2mm	88.023/2	88.024/2	88.025/2	
	3mm	88.023/3	88.024/3	88.025/3	
	4mm	88.023/4	88.024/4	88.025/4	
	5mm				88.025/5

To see the DESSLoc® range please refer to page 53.

Scan Bodys					
	for desktop scanning * use indicated screw (not included)		50.023 * 19.015	50.024 * 19.016	50.025 * 19.017
	ON UNIABUTMENT®				
	20°	50.036 * 19.018			
	45°	50.037 * 19.018			
for intra-oral scanning			52.023	52.024	52.025
ON UNIABUTMENT®					
20°		52.036			
45°		52.037			
on Ti base		●	53.100-P10	53.100-P10	53.100-P10
Pack of 10		○	53.200-P10	53.200-P10	53.200-P10

Pre-milled blank Abutments

YELLOW

AQUA

LILAC

3.0

3.5-4.0

4.5-5.0

	Ti	ø10mm	61.023	61.024	61.025
		ø14mm		62.024	62.025
	CoCr	ø14mm		71.024	71.025

UniAbutment®

	20°	0,5mm		40.013/0	40.015/0
		1,0mm		40.013/1	40.015/1
		2,0mm		40.013/2	40.015/2
		4,0mm		40.013/4	40.015/4
		6,0mm		40.013/6	40.015/6
	45°	0,5mm		40.014/0	40.016/0
		1,0mm		40.014/1	40.016/1
		2,0mm		40.014/2	40.016/2
		4,0mm		40.014/4	40.016/4
		6,0mm		40.014/6	40.016/6

Multi-unit abutments

MULTI-UNIT PLATFORM

4.8

	Straight	1,5mm	●	40.024/1	40.025/1
		2,5mm	●	40.024/2	40.025/2
		3,5mm	●	40.024/3	40.025/3
	Angled 17°	3,5mm	●	41.724/3	41.725/3
		4,5mm	●	41.724/4	41.725/4
	30°	4,5mm	●	43.024/4	43.025/4
		5,5mm	●	43.024/5	43.025/5

To see the components on top of the DESS Multiunits please refer to page 48.



INTERNAL HEX MI

compatible with
MIS® SEVEN

All abutments include final screw, except all plastic castables in packs of 10 and scan bodys for desktop scanning.

- non-engaging
- engaging

		NP	SP	WP	Multiunit	
		YELLOW	PURPLE	GREEN	4.8	
		3.3	3.75-4.2	5.0-6.0		
	Single	19.056	19.081	19.081	19.110	
	Single				19.085	
	Pack of 10		19.081-P10	19.081-P10		

Titanium Analog - For 3D Printed Models

	Single	14.056	14.095	14.096	14.007	
	Single		(for MIS MUA only)		14.085	
	Pack of 10		14.095-P10	14.096-P10	14.007 P10	

Our analogs are suitable for both 3D and traditional models.

Straight Ti base

New

	0,3mm				15.096	15.085	
	0,8mm		15.056	15.095			
	1,5mm		15.256	15.295	15.296		
	3,0mm		15.356	15.395	15.396		
	0,3mm				16.096		
	0,8mm		16.056	16.095			
	1,5mm		16.256	16.295	16.296		
	3,0mm		16.356	16.395	16.396		

Ti base Castables

	Pack of 5		32.001-P5	32.001-P5	32.001-P5
			33.001-P5	33.001-P5	33.001-P5

 For DESS® Multi-Unit

 For MIS® Multi-Unit

AngleBase®



Plastic sleeves for AngleBase®

Pack of 5	0°	33.100-P5	33.100-P5
Pack of 5	10°	33.101-P5	33.101-P5
Pack of 5	20°	33.102-P5	33.102-P5

Torx® Screw for AngleBase®



Torx® ball Screwdriver	20mm	24mm	30mm	35mm
	DT20TB06	DT24TB06	DT30TB06	DT35TB06

to use with AngleBase® & ELLIPTIBase®

CoCr Abutment



Includes screw and wax up sleeve

	30.056	30.095	30.096
	31.056	31.095	31.096

All Plastic Castable Abutments



Single		28.056	28.095	28.096
Pack of 10			11.017-P10	11.018-
Single		29.056	29.095 P10	29.096
Pack of 10			12.017-P10	12.018-P10

Impression Transfer	NP	SP	WP
	YELLOW	PURPLE	GREEN
	3.3	3.75-4.2	5.0-6.0
	17.056	17.017	17.018



Impression transfer packaged with open and closed tray screw and transfer.

Titanium Temporary Abutments

●	24.095	24.096
●	25.095	25.096



Abutments



Straight	0°	13.095	13.096
	15°	22.095	22.096
	25°	23.095	23.096

Healing Abutments



3,0mm	20.017/3	20.018/3
4,0mm	20.017/4	20.018/4
5,0mm	20.017/5	20.018/5

DESSLoc®



1,0mm	88.017/1	88.018/1
2,5mm	88.017/2	88.018/2
3,5mm	88.017/3	88.018/3
4,5mm	88.017/4	88.018/4
5,5mm	88.017/5	

To see the DESSLoc® range please refer to page 53.

Scan Bodys



Scan Bodys	NP	SP	WP	Multiunit
	YELLOW	PURPLE	GREEN	4.8
	3.3	3.75-4.2	5.0-6.0	
for desktop scanning	50.056 * 19.056	50.017 * 19.081	50.018 * 19.081	
* use indicated screw (not incuded)				
for intra-oral scanning	52.056	52.017	52.018	52.085 MIS
on Ti base				
Pack of 10	● 53.100-P10	53.100-P10	53.100-P10	
	● 53.200-P10	53.200-P10	53.200-P10	

Pre-milled blank Abutments



Ti	ø10mm	61.056	61.095	61.096
CoCr	ø14mm		71.095	71.096

To see all the Pre-milled blank available please refer to page 138.

Multi-unit abutments



				MULTI-UNIT PLATFORM 4.8	
Straight	1,0mm	●		40.019/1	40.020/1
	2,0mm	●		40.019/2	40.020/2
	3,0mm	●		40.019/3	40.020/3
	4,0mm	●		40.019/4	
	5,0mm	●		40.019/5	
Angled	17°	2,5mm	●	41.719/2	41.720/2
		3,5mm	●	41.719/3	41.720/3
	30°	3,5mm	●	43.019/3	43.020/3
		4,5mm	●	43.019/4	43.020/4

MIS For MIS® Multi-Unit

To see the components on top of the DESS Multiunits please refer to page 48.



CONIC MC

compatible with
MIS® C1 internal

All abutments include final screw

-  non-engaging
-  engaging

		NP	SP	WP	Multi-unit
Final Screws - .050/1.27Hex		YELLOW	PURPLE	GREEN	4.8
	Single	19.056	19.081	19.081	19.085 

Titanium Analog - For 3D Printed Models					
	Single	14.081	14.082	14.083	14.085 

Our analogs are suitable for both 3D and traditional models.

Straight Ti base		GH			
	0,3mm				15.085 
	1,0mm		15.081	15.082	15.083
	1,0mm		16.081	16.082	16.083
		Ti base Castables			
	Pack of 5		32.001-P5	32.001-P5	32.001-P5
			33.001-P5	33.001-P5	33.001-P5

Scan Bodys					
	for intra-oral scanning		52.081	52.082	52.083
					52.085 
	on Ti base		53.100-P10	53.100-P10	53.100-P10
	Pack of 10		53.200-P10	53.200-P10	53.200-P10

Pre-milled blank Abutments					
	Ti	ø10mm	61.081	61.082	61.083

New

 For MIS® Multi-Unit



CONIC MV

compatible with
MIS® V3

All abutments include final screw

- non-engaging
- engaging

		NP	SP	Multi-unit	 For MIS® Multi-Unit
Final Screws - .050/1.27Hex		blue	purple	4.8	
	Single	19.056	19.081	19.085	

Titanium Analog - For 3D Printed Models					
	Single	14.084	14.082	14.085	

Our analogs are suitable for both 3D and traditional models.

Straight Ti base		GH			
	0,3mm	●			15.085 
	1,0mm	●	15.084	15.082	
	1,0mm	●	16.084	16.082	
			Ti base Castables		
	Pack of 5	●	32.001-P5	32.001-P5	32.007-P5
		●	33.001-P5	33.001-P5	

Scan Bodys					
	for intra-oral scanning		52.084	52.082	52.085 
	on Ti base				
	Pack of 10	●	53.100-P10	53.100-P10	
		●	53.200-P10	53.200-P10	

Pre-milled blank Abutments					
	Ti	ø10mm	61.084	61.082	

 For MIS® Multi-Unit

New

30



CONIC ANYR

compatible with
MEGAGEN AnyRidge®

All abutments include final screw, except scan bodys for desktop scanning.

Final Screws - .048/1.20 Hex



Single	19.057
--------	--------

Titanium Analog - For 3D Printed Models



Single	14.057
--------	--------

Our analogs are suitable for both 3D and traditional models.

Straight Ti base



GH		
2,0mm	●	15.057
3,0mm	●	15.357
2,0mm	●	16.057
3,0mm	●	16.357

Ti base Castables



Pack of 5	●	32.001-P5
	●	33.001-P5

- non-engaging
- engaging

AngleBase®



●	36.057
●	37.057

Plastic sleeves for AngleBase®



Pack of 5	0°	33.100-P5
Pack of 5	10°	33.101-P5
Pack of 5	20°	33.102-P5

Torx® Screw for AngleBase®



(replacement screw)

19.460

Torx® ball Screwdriver	20mm	24mm	30mm	35mm
	DT20TB06	DT24TB06	DT30TB06	DT35TB06

to use with AngleBase® & ELLIPTIBase®

Scan Bodys



for desktop scanning	50.057
* use indicated screw (not incuded)	* 19.057
for intra-oral scanning	52.057
on Ti base	
Pack of 10	 53.100-P10
	 53.200-P10

Pre-milled blank Abutments



Ti	ø10mm	61.057
CoCr	ø14mm	71.057

To see all the Pre-milled blank available please refer to page 138.



CONIC OSS

compatible with
HIOSSEN® ET

-  non-engaging
-  engaging

All abutments include final screw.

			Mini 3.5-4.0 Yellow	Standard 4.5-5.0 Green
	Single	DLC coated	19.072	19.073
	Single		19.272	19.273

Titanium Analog - For 3D Printed Models

	Single		14.072	14.073
Our analogs are suitable for both 3D and traditional models.				

Straight Ti base

GH				
	1,0mm		15.072	15.073
	3,0mm		15.372	15.373
	1,0mm		16.072	16.073
	3,0mm		16.372	16.373

Ti base Castables

	Pack of 5		32.001-P5	32.001-P5
			33.001-P5	33.001-P5

C-Base® - For Cerec

	1,0mm		58.072	58.073
	3,0mm		58.072/3	58.073/3
	1,0mm		59.072	59.073
	3,0mm		59.072/3	59.073/3

Includes screw, same as standard screw.

* Either the .048/1.20 hex or 1.22 hex driver will work with the Hiossen screws.

		Mini 3.5-4.0 Yellow	Standard 4.5-5.0 Green
		36.072	36.073
		37.072	37.073

Plastic sleeves for AngleBase®

	Pack of 5	0°	33.100-P5	33.100-P5
	Pack of 5	10°	33.101-P5	33.101-P5
	Pack of 5	20°	33.102-P5	33.102-P5

Torx® Screw for AngleBase®

	(replacement screw)	19.444	19.461
---	---------------------	--------	--------

Torx® ball Screwdriver	20mm	24mm	30mm	35mm
	DT20TB06	DT24TB06	DT30TB06	DT35TB06

to use with AngleBase® & ELLIPTIBase®

CoCr Abutment

	Includes screw and wax up sleeve		30.072	30.073
			31.072	31.073

Titanium Temporary Abutments

		24.072	24.073
		25.072	25.073

New

New

New

Scan Bodys



		Mini 3.5-4.0 Yellow	Standard 4.5-5.0 Green
for desktop scanning		50.072	50.073
* use indicated screw (not included)		* 19.072	* 19.073
for intra-oral scanning		52.072	52.073
on Ti base	Pack of 10	53.100-P10	53.100-P10
		53.200-P10	53.200-P10

Pre-milled blank Abutments



Ti	ø10mm	61.072	61.073
	ø14mm	62.072	62.073
CoCr	ø14mm	71.072	71.073

Multi-unit abutments



		MULTI-UNIT PLATFORM 4.8	
Straight	1,5mm	40.072/1	40.073/1
	2,5mm	40.072/2	40.073/2
	3,5mm	40.072/3	40.073/3



Angled 17°	2,5mm	41.072/2	41.073/2
	3,5mm	41.072/3	41.073/3



30°	3,5mm	43.072/3	43.073/3
	4,5mm	43.072/4	43.073/4

To see the components on top of the DESS Multiunits please refer to page 48.



INTERNAL CAM

compatible with
Camlog®

All abutments include final screw, except castables in packs of 10 and scan abutments for desktop scanning.

- non-engaging
- engaging

		SILVER	YELLOW	PINK	BLUE
Final Screws - .050/1.27Hex		3.3	3.8	4.3	5.0
	Single	19.044	19.044	19.044	19.045

Titanium Analog - For 3D Printed Models

	Single	14.064	14.065	14.066	14.067
Our analogs are suitable for both 3D and traditional models.					

Straight Ti base

		GH			
	0,3mm	●		15.065	15.066
	0,5mm	●	15.064		15.067
	1,5mm	●		15.265	15.266
	3,0mm	●		15.365	15.366
	0,3mm	●		16.065	16.066
	0,5mm	●	16.064		16.067
	1,5mm	●		16.265	16.266
	3,0mm	●		16.365	16.366



Ti base Castables

Pack of 5	●	32.001-P5	32.001-P5	32.001-P5	32.001-P5
	●	33.001-P5	33.001-P5	33.001-P5	33.001-P5

CoCr Abutment

		30.065	30.066
		31.065	31.066
Includes screw and wax up sleeve			

AngleBase®



YELLOW PINK

3.8 4.3

-
-

36.065 36.066
37.065 37.066

Plastic sleeves for AngleBase®



Pack of 5	0°	33.100-P5	33.100-P5
Pack of 5	10°	33.101-P5	33.101-P5
Pack of 5	20°	33.102-P5	33.102-P5

Torx® Screw for AngleBase®



(replacement screw)

19.444 19.444

Torx® ball Screwdriver	20mm	24mm	30mm	35mm
	DT20TB06	DT24TB06	DT30TB06	DT35TB06

to use with AngleBase® & ELLIPTIBase®

All Plastic Castable Abutments



-

28.065 28.066

Includes final screw

Transfer



17.065 17.066

Impression transfer packaged with open and closed tray screw and transfer.

Titanium Temporary Abutments

SILVER

YELLOW

PINK

BLUE

3.3

3.8

4.3

5.0



24.065	24.066
25.065	25.066

Healing Abutments



3,0mm	20.065/3	20.066/3
4,0mm	20.065/4	20.066/4
5,0mm	20.065/5	20.066/5

DESSLoc®



1,0mm	88.065/1	88.066/1
2,0mm	88.065/2	88.066/2
3,0mm	88.065/3	88.066/3
4,0mm	88.065/4	88.066/4

To see the DESSLoc® range please refer to page 53.

Scan Bodys



for desktop scanning	50.064	50.065	50.066	50.067
* use indicated screw (not included)	* 19.044	* 19.044	* 19.044	* 19.045
for intra-oral scanning	52.064	52.065	52.066	52.067
on Ti base	53.100-P10	53.100-P10	53.100-P10	53.100-P10
Pack of 10	53.200-P10	53.200-P10	53.200-P10	53.200-P10

Pre-milled blank Abutments



Ti	ø10mm	61.064	61.065	61.066	61.067
	ø14mm		62.065	62.066	62.067
CoCr	ø14mm	71.064	71.065	71.066	71.067



INTERNAL CON

compatible with
Conelog®

All abutments include final screw

Final Screws - .050/1.27Hex		SILVER	YELLOW	BLUE
		3.3	3.8/4.3	5.0
	Single	19.074	19.074	19.076

Titanium Analog - For 3D Printed Models		SILVER	YELLOW	BLUE
	Single	14.074	14.075	14.076
Our analogs are suitable for both 3D and traditional models.				

Straight Ti base		GH			
	1,4mm	●	15.074	15.075	15.076
	1,4mm	●	16.074	16.075	16.076
Ti base Castables					
	Pack of 5	●	32.001-P5	32.001-P5	32.001-P5
		●	P5	33.001-P5	33.001-P5

C-Base®					
	1,0mm	●	58.074	58.075	58.076
	2,1mm	●	58.074/2	58.075/2	58.076/2
	1,0mm	●	59.074	59.075	59.076
	2,1mm	●	59.074/2	59.075/2	59.076/2

Replacement Screw for C-Base®				
		19.052	19.052	19.053

- non-engaging
- engaging



Torx® Screw for AngleBase®					
	(replacement screw)		19.474	19.474	19.476
Torx® ball Screwdriver		20mm	24mm	30mm	35mm
		DT20TB06	DT24TB06	DT30TB06	DT35TB06
to use with AngleBase® & ELLIPTIBase®					

Scan Bodys					
	for intra-oral scanning		52.074	52.075	52.076
	on Ti base	Pack of 10	●	53.100-P10	53.100-P10
			●	53.200-P10	53.200-P10

Pre-milled blank Abutments					
	Ti	ø10mm	61.074	61.075	61.076
		ø14mm	62.074	62.075	62.076
	CoCr	ø14mm	71.074	71.075	71.076



INTERNAL ANK

compatible with
Dentsply Ankylos[®] C/X

All abutments include final screw, except castables in packs of 10 and scan abutments for desktop scanning.

- non-engaging
- engaging

Final Screw - 1.00 Hex

	height:		
	Single	1,5mm	19.048
	Pack of 10		19.048-P10
	Single	2,5mm	19.050

Titanium Analog - For 3D Printed Models

	Single	Ø 3,5	14.048
	Pack of 10		14.048-P10
	Single	Ø 4,5	14.049
	Pack of 10		14.049-P10

Our analogs are suitable for both 3D and traditional models.

Straight Ti base

	GH		
	1,5mm	●	15.050
	2,5mm	●	15.250
	1,5mm	○	16.050
	2,5mm	○	16.250

Ti base Castables



- 32.001-P5
- 33.001-P5

C-Base[®] - For Cerec

	1,4mm	●	58.050
	2,5mm	●	58.050/2
	1,4mm	○	59.050
	2,5mm	○	59.050/2

Replacement Screw for C-Base[®]

 (replacement screw)	19.046
---	--------

CoCr Abutment

	height:		
	1,5mm	●	30.050
	2,5mm	●	30.250
	1,5mm	○	31.050
	2,5mm	○	31.250

Includes screw and wax up sleeve

All Plastic Castable Abutments

	1,5mm	●	28.048
	2,5mm	●	28.250

Impression Transfer

	1,5mm	○	17.049
	2,5mm	○	17.250

Impression transfer packaged with open and closed tray screw and transfer.

Titanium Temporary Abutments

	1,5mm	●	24.048
	2,5mm	●	24.250
	1,5mm	○	25.048
	2,5mm	○	25.250

New

New

Abutments



		height:		
15°	Straight	1,5mm		13.047
		2,5mm		13.250
	Angled	1,5mm		22.048
		2,5mm		22.250
25°	Angled	1,5mm		23.048
		2,5mm		23.250

Healing Abutments



		height:	
		1,5mm	20.046/1
		2,5mm	20.046/3

DESSLoc®



		height:	
		2,0mm	88.046/2
		3,0mm	88.046/3
		4,0mm	88.046/4

To see the DESSLoc® range please refer to page 53.

Scan Bodys



for desktop scanning		50.049
* use indicated screw (not included)		* 19.048
for intra-oral scanning		52.049
on Ti base	Pack of 10	53.100-P10
		53.200-P10

Pre-milled blank Abutments



Ti	ø10mm	61.250
CoCr	ø14mm	71.250



CONIC DENT

compatible with
Dentium Superline™ & Implantium®

All abutments include final screw, except castables in packs of 10 and scan abutments for desktop scanning.

- non-engaging
- ⬡ engaging

Final Screw - .050/1.27 hex

RP / YELLOW / GREEN / BLUE / RED
3.8 3.6 4.0 4.5 5.0

	Single	19.090
	ON ABUTMENT	
	Single	19.066

Titanium Analog - For 3D Printed Models



Our analogs are suitable for both 3D and traditional models.

Straight Ti base

GH		
0,9mm	●	15.090
1,5mm	●	15.190
3,0mm	●	15.390
0,9mm	⬡	16.090
1,5mm	⬡	16.190
3,0mm	⬡	16.390

Ti base Castables

	Pack of 5	●	32.001-P5
		⬡	33.001-P5

AngleBase®



RP / YELLOW / GREEN / BLUE / RED
3.8 3.6 4.0 4.5 5.0

●	36.090
⬡	37.090

New

Plastic sleeves for AngleBase®



Pack of 5	0°	33.100-P5
Pack of 5	10°	33.101-P5
Pack of 5	20°	33.102-P5

Torx® Screw for AngleBase®



(replacement screw)

19.405

New

Torx® ball Screwdriver 20mm 24mm 30mm 35mm



DT20TB06 DT24TB06 DT30TB06 DT35TB06

to use with AngleBase® & ELLIPTIBase®

CoCr Abutment



●	30.090
⬡	31.090

New

Includes screw and wax up sleeve

Impression Transfer



Impression transfer packaged with open and closed tray screw and transfer.

RP / YELLOW / GREEN / BLUE / RED
3.8 3.6 4.0 4.5 5.0

17.090

Scan Bodys



for desktop scanning
* use indicated screw (not incuded) 50.090
* 19.090

for intra-oral scanning 52.090

on Ti base 53.100-P10
Pack of 10 53.200-P10

Pre-milled blank Abutments



Ti	ø10mm	61.090
	ø14mm	62.090
CoCr	ø14mm	71.090

To see all the Pre-milled blank available please refer to page 138.

Multi-unit abutments



New

Straight	1,2mm	●	40.090/1	
	2,5mm	●	40.090/2	
	3,5mm	●	40.090/3	
Angled	17°	2,5mm	◆	41.790/2
		3,5mm	◆	41.790/3
	30°	2,5mm	◆	43.090/2
		3,5mm	◆	43.090/3

To see the components on top of the DESS Multiunits please refer to page 48.



INTERNAL HEX FD

compatible with
Xive®

All abutments include final screw, except castables in packs of 10 and scan bodys for desktop scanning.

- non-engaging
- ⬡ engaging

		SILVER	YELLOW	BLUE	RED
		3.4	3.8	4.5	5.5
	Single	19.024	19.024	19.024	19.024
	Pack of 10	19.024-P10	19.024-P10	19.024-P10	19.024-P10

Titanium Analog - For 3D Printed Models

	Single	14.038	14.039	14.040	14.071
	Pack of 10		14.039-P10		
Our analogs are suitable for both 3D and traditional models.					

Straight Ti base

	0,3mm	●	15.038	15.039	15.040	15.071
	1,5mm	●	15.238	15.239	15.240	15.271
	3,0mm	●	15.338	15.339	15.340	15.371
	0,3mm	⬡	16.038	16.039	16.040	16.071
	1,5mm	⬡	16.238	16.239	16.240	16.271
	3,0mm	⬡	16.338	16.339	16.340	16.371

Ti base Castables

	Pack of 5	●	32.001-P5	32.001-P5	32.001-P5	32.001-P5
		⬡	33.001-P5	33.001-P5	33.001-P5	33.001-P5

			SILVER	YELLOW	BLUE	RED
			3.4	3.8	4.5	5.5
	0,3mm	●	58.038	58.039		
	0,6mm	●			58.040	58.071
	1,5mm	●	58.038/2	58.039/2	58.040/2	58.071/2
	3,0mm	●	58.038/3	58.039/3	58.040/3	58.071/3
	0,3mm	⬡	59.038	59.039		
	0,6mm	⬡			59.040	59.071
	1,5mm	⬡	59.038/2	59.039/2	59.040/2	59.071/2
	3,0mm	⬡	59.038/3	59.039/3	59.040/3	59.071/3
	Includes screw, same as standard screw.		Please refer to page 156 for additional information.			

New

AngleBase®

	●	36.038	36.039	36.040
	⬡	37.038	37.039	37.040

Plastic sleeves for AngleBase®

	Pack of 5	0°	33.100-P5	33.100-P5
	Pack of 5	10°	33.101-P5	33.101-P5
	Pack of 5	20°	33.102-P5	33.102-P5

Torx® Screw for AngleBase®

	(replacement screw)	19.438	19.438	19.438	
Torx® ball Screwdriver		20mm	24mm	30mm	35mm
		DT20TB06	DT24TB06	DT30TB06	DT35TB06


to use with AngleBase® & ELLIPTIBase®

		SILVER	YELLOW	BLUE	RED
		3.4	3.8	4.5	5.5
		●	30.038	30.039	30.040
		●	31.038	31.039	31.040
Includes screw and wax up sleeve					

All Plastic Castable Abutments

	Single	●	28.038	28.039	28.040
	Pack of 10		11.039-P10		
	Single	●	29.038	29.039	29.040
	Pack of 10		12.039-P10		

Impression Transfer

			17.038	17.039	17.040
	Impression transfer packaged with open and closed tray screw and transfer.				

Titanium Temporary Abutments

		●	24.038	24.039	24.040
		●	25.038	25.039	25.040

Abutments

	Straight		13.038	13.039	13.040
	Angled	15°	22.038	22.039	22.040
		25°		23.039	

Healing Abutments

	3,0mm		20.036/3	20.037/3	20.038/3
	4,0mm		20.036/4	20.037/4	20.038/4
	5,0mm		20.036/5	20.037/5	20.038/5

		SILVER	YELLOW	BLUE	RED
		3.4	3.8	4.5	5.5
	1,0mm	88.036/1	88.037/1		
	2,0mm	88.036/2	88.037/2		
	3,0mm	88.036/3	88.037/3		
	4,0mm	88.036/4	88.037/4		
To see the DESSLoc® range please refer to page 53.					

Scan Bodies

	for desktop scanning	50.038	50.039	50.040	50.071
	* use indicated screw (not included)	* 19.024	* 19.024	* 19.024	* 19.024
	for intra-oral scanning	52.038	52.039	52.040	52.071
	on Ti base				
Pack of 10	●	53.100-P10	53.100-P10	53.100-P10	53.100-P10
	●	53.200-P10	53.200-P10	53.200-P10	53.200-P10

Pre-milled blank Abutments

	Ti	ø10mm	61.038	61.039	61.040	61.071
		ø14mm		62.039	62.040	62.071
	CoCr	ø14mm	71.038	71.039	71.040	71.071



EXTERNAL HEX UNIVERSAL

compatible with
Nobel Brånemark®

All abutments include final screw, except castables in packs of 10 and scan bodys for desktop scanning.

- non-engaging
- ⬡ engaging

		NP	RP	WP
Final Screws - Unigrip		3.5	4.1	5.1
	Single	19.027	19.028	19.029
	Pack of 10	19.027-P10	19.028-P10	19.029-P10
	DLC coated	19.227	19.228	19.229

Titanium Analog - For 3D Printed Models		NP	RP	WP
	Single	14.001	14.002	14.003
	Pack of 10		14.002-P10	
Our analogs are suitable for both 3D and traditional models.				

Straight Ti base		GH			
	0,3mm	●	15.001	15.002	15.003
	1,5mm	●	15.201	15.202	15.203
	3,0mm	●	15.301	15.302	15.303
	0,3mm	⬡	16.001	16.002	16.003
	1,5mm	⬡	16.201	16.202	16.203
	3,0mm	⬡	16.301	16.302	16.303
Ti base Castables					
	Pack of 5	●	32.001-P5	32.001-P5	32.003-P5
		⬡	33.001-P5	33.001-P5	33.003-P5

C-Base®					
	1,5mm	●	58.001/2	58.002/2	
	3,0mm	●	58.001/3	58.002/3	
	0,6mm	⬡	59.001	59.002	
	1,5mm	⬡	59.001/2	59.002/2	
	3,0mm	⬡	59.001/3	59.002/3	
Includes screw, same as standard screw.					

New

AngleBase®			NP	RP	WP
			3.5	4.1	5.1
	0,3mm	●	36.001	36.002	
	1,5mm	●	36.001/2	36.002/2	
	2,5mm	●	36.001/3	36.002/3	
	0,3mm	⬡	37.001	37.002	
	1,5mm	⬡	37.001/2	37.002/2	
	2,5mm	⬡	37.001/3	37.002/3	

Plastic sleeves for AngleBase®					
	Pack of 5	0°	33.100-P5	33.100-P5	
	Pack of 5	10°	33.101-P5	33.101-P5	
	Pack of 5	20°	33.102-P5	33.102-P5	

Torx® Screw for AngleBase®					
	(replacement screw)		19.401	19.402	
Torx® ball Screwdriver					
	20mm	DT20TB06	DT24TB06	DT30TB06	DT35TB06
to use with AngleBase® & ELLIPTIBase®					

CoCr Abutment					
		●	30.001	30.002	30.003
		⬡	31.001	31.002	31.003
Includes screw and wax up sleeve					

			NP	RP	WP
All Plastic Castable Abutments			3.5	4.1	5.1
	Single	●	28.001	28.002	28.003
	Pack of 10	●		11.002-P10	
	Single	●	29.001	29.002	29.003
	Pack of 10	●		12.002-P10	

Impression Transfer

			17.041	17.042	17.043
Impression transfer packaged with open and closed tray screw and transfer.					

Titanium Temporary Abutments

	●	24.001	24.002	24.003
	●	25.001	25.002	25.003


Abutments

	Straight		13.001	13.002	13.003
	Angled 15°		22.001	22.002	22.003
	Angled 25°			23.002	

Healing Abutments

	3,0mm	20.001/3	20.002/3	20.003/3
	4,0mm	20.001/4	20.002/4	20.003/4
	5,0mm	20.001/5	20.002/5	20.003/5

DESSLoc®

	1,0mm	88.001/1	88.002/1	
	2,0mm	88.001/2	88.002/2	
	3,0mm	88.001/3	88.002/3	
	4,0mm		88.002/4	
	5,0mm		88.002/5	

To see the DESSLoc® range please refer to page 53.

Scan Bodys

			NP	RP	WP
			3.5	4.1	5.1
		for desktop scanning * use indicated screw- sold seperately	50.001 * 19.027	50.002 * 19.028	50.003 * 19.029
		for intra-oral scanning	52.001	52.002	52.003
on Ti base					
		●	53.100-P10	53.100-P10	53.100-P10
Pack of 10		●	53.200-P10	53.200-P10	53.200-P10

Pre-milled blank Abutments

	Ti	ø10mm	61.001	61.002	61.003
	CoCr	ø14mm	71.001	71.002	71.003
To see all the Pre-milled blank available please refer to page 138.					

Multi-unit abutments

		MULTI-UNIT PLATFORM			
		4.8		6.5	
	Straight	1,0mm	●	40.001/1	40.002/1 40.003/1
		2,0mm	●	40.001/2	40.002/2 40.003/2
		3,0mm	●	40.001/3	40.002/3 40.003/3
		4,0mm	●		40.002/4
		5,0mm	●		40.002/5
	Angled 17°	2,0mm	●	41.701/2	41.702/2
		3,0mm	●	41.701/3	41.702/3
		4,0mm	●		41.702/4
		30° 4,0mm	●		43.002/4
		5,0mm	●		43.002/5

All Multi-units include DLC SCREW.

To see the components on top of the DESS Multiunits please refer to page 48.



EXTERNAL HEX USA

compatible with
3i Osseotite® External Hex

All abutments include final screw, except all plastic castables in packs of 10 and scan bodys for desktop scanning.

- non-engaging
- ⬡ engaging

Final Screws - 048/1.20 Hex		PURPLE	BLUE	YELLOW
		3.4	4.1	5.0



Single	19.034	19.034	19.034
Pack of 10	19.034-P10	19.034-P10	19.034-P10

Titanium Analog - For 3D Printed Models



Single	14.011	14.012	14.013
Pack of 10		14.012-P10	

Our analogs are suitable for both 3D and traditional models.

Straight Ti base



GH				
0,3mm	●	15.011	15.012	15.013
1,5mm	●	15.211	15.212	15.213
3,0mm	●	15.311	15.312	15.313
0,3mm	⬡	16.011	16.012	16.013
1,5mm	⬡	16.211	16.212	16.213
3,0mm	⬡	16.311	16.312	16.313

Ti base Castables



Pack of 5	●	32.001-P5	32.001-P5	32.001-P5
	⬡	33.001-P5	33.001-P5	33.001-P5

C-Base®



1,5mm	●	58.011/2	58.012/2	58.013/2
3,0mm	●	58.011/3	58.012/3	58.013/3
0,6mm	⬡	59.011	59.012	59.013
1,5mm	⬡	59.011/2	59.012/2	59.013/2
3,0mm	⬡	59.011/3	59.012/3	59.013/3

Includes screw, same as standard screw.

Please refer to page 156 for additional information.

New

AngleBase®



		PURPLE	BLUE	YELLOW
		3.4	4.1	5.0

0,3mm	●	36.011	36.012	36.013
1,5mm	●	36.011/2	36.012/2	36.013/2
2,5mm	●	36.011/3	36.012/3	36.013/3
0,3mm	⬡	37.011	37.012	37.013
1,5mm	⬡	37.011/2	37.012/2	37.013/2
2,5mm	⬡	37.011/3	37.012/3	37.013/3

Plastic sleeves for AngleBase®



Pack of 5	0°	33.100-P5	33.100-P5	33.100-P5
Pack of 5	10°	33.101-P5	33.101-P5	33.101-P5
Pack of 5	20°	33.102-P5	33.102-P5	33.102-P5

Torx® Screw for AngleBase®



(replacement screw)

19.402	19.402	19.402
--------	--------	--------

Torx® ball Screwdriver	20mm	24mm	30mm	35mm
	DT20TB06	DT24TB06	DT30TB06	DT35TB06



CoCr Abutment



Includes screw and wax up sleeve

●	30.011	30.012	30.013
⬡	31.011	31.012	31.013

to use with AngleBase® & ELLIPTIBase®

All Plastic Castable Abutments

PURPLE 3.4BLUE 4.1YELLOW 5.0

	Single		●	28.011	28.012	28.013
	Pack of 10			11.002-P10		
	Single		◆	29.011	29.012	29.013
	Pack of 10			12.002-P10		

Impression Transfer

				17.011	17.012	17.013
	Impression Impression transfer packaged with open and closed tray screw and transfer.					

Titanium Temporary Abutments

			●	24.011	24.012	24.013
			◆	25.011	25.012	25.013

Abutments

	Straight			13.011	13.012	13.013
	Angled	15°		22.011	22.012	22.013
		25°			23.012	

Healing Abutments

	3,0mm	20.011/3	20.012/3	20.013/3
	4,0mm	20.011/4	20.012/4	20.013/4
	5,0mm	20.011/5	20.012/5	20.013/5

DESSLoc®

	1,0mm			88.011/1	88.002/1	
	2,0mm			88.011/2	88.002/2	
	3,0mm			88.011/3	88.002/3	
	4,0mm			88.011/4	88.002/4	
	5,0mm				88.002/5	

To see the DESSLoc® range please refer to page 53.

Scan Bodys

	for desktop scanning * use indicated screw- sold seperately			50.011 * 19.034	50.012 * 19.034	50.013 * 19.034
	for intra-oral scanning			52.011	52.012	52.013
	on Ti base	Pack of 10	●	53.100-P10	53.100-P10	53.100-P10
			◆	53.200-P10	53.200-P10	53.200-P10

Pre-milled blank Abutments

	Ti	ø10mm		61.011	61.012	61.013
	CoCr	ø14mm		71.011	71.012	71.013
To see all the Pre-milled blank available please refer to page 138.						

Multi-unit abutments

			MULTI-UNIT PLATFORM 4.8			
	Straight	1,0mm	●	40.009/1	40.017/1	
		2,0mm	●	40.009/2	40.017/2	
		3,0mm	●	40.009/3	40.017/3	
		4,0mm	●		40.017/4	
		5,0mm	●		40.017/5	
	Angled 17°	2,0mm	◆		41.712/2	
		3,0mm	◆		41.712/3	
		4,0mm	◆		41.712/4	
	30°	3,0mm	◆		43.012/3	
		4,0mm	◆		43.012/4	
		5,0mm	◆		43.012/5	

To see the components on top of the DESS Multiunits please refer to page 48.



EXTERNAL HEX BH

compatible with
BioHorizons® External Hex

All abutments include final screw, except castables in packs of 10 and scan bodys for desktop scanning.

- non-engaging
- engaging

		YELLOW 3.7	GREEN 4.2	BLUE 5.2
	Single	19.013	19.014	19.014
	Pack of 10		19.014-P10	19.014-P10

Titanium Analog - For 3D Printed Models

	Single	14.020	14.021	14.022
	Pack of 10		14.021-P10	14.022-P10

Our analogs are suitable for both 3D and traditional models.

Straight Ti base		GH			
	0,3mm				
	0,3mm				

Ti base Castables

All Plastic Castable Abutments

	Single			
	Pack of 10			
	Single			
	Pack of 10			

Impression Transfer

Impression transfer packaged with open and closed tray screw and transfer.

Titanium Temporary Abutments



Abutments



Healing Abutments



Scan Bodys



Pre-milled blank Abutments



YELLOW 3.7 GREEN 4.2 BLUE 5.2

		24.020	24.021	24.022
		25.020	25.021	25.022

Straight		13.020	13.021	13.022
Angled 15°		22.020	22.021	

3mm		20.020/3	20.021/3	20.022/3
4mm		20.020/4	20.021/4	20.022/4
5mm		20.020/5	20.021/5	20.022/5

for desktop scanning * use indicated screw- sold seperately		50.020 * 19.013	50.021 * 19.014	50.022 * 19.014
--	--	--------------------	--------------------	--------------------

for intra-oral scanning		52.020	52.021	52.022
-------------------------	--	--------	--------	--------

on Ti base				
Pack of 10		53.100-P10	53.100-P10	53.100-P10
		53.200-P10	53.200-P10	53.200-P10

Ti	ø10mm	61.020	61.021	61.022
		62.020	62.021	62.022
CoCr	ø14mm	71.020	71.021	71.022



MULTIUNIT CONNECTION

All of these "RP" components will work on top of the following multiunit abutments - Nobel, BioHorizons, 3i, Megagen, Hiossen, Neodent, BSB, Astra EV Multibase, and DESS Multiunits. All are 4.8mmD.

The "WP" components are only for Branemark WP multiunits only.

Final Screws			RP	WP
	Nobel Unigrip	Single	19.006	19.031
		Pack of 10	19.006-P10	19.031-P10
		with DLC-coating	19.206	19.231
	.048/1.2 Hex	Single	19.055	
	.050/1.27 Hex	Single	19.110	

All of these screws will work with all the "RP" components, except the AngleBase.

Analog			
	Single	14.007	14.047
	Pack of 10	14.007-P10	14.047-P10

Our analogs are for both 3D and traditional models.

Straight Ti base		Shaft height	
	6,0mm	15.007	15.054
	8,0mm	15.007/8	

Includes screw 19.006 for 15.007 and 15.007/8 and screw 19.031 for 15.054

Ti base Castables			
	Pack of 5	32.007-P5	32.054-P5
	This is for usage with part 15.007 only		

AngleBase®		RP	WP
		36.007	36.054
	Includes screw 19.306 (for RP) and 19.331 (for WP)		
Torx® Screw for AngleBase®		19.306	19.331
	Torx® ball Screwdriver	20mm	24mm
		DT20TB06	DT24TB06
		30mm	35mm
		DT30TB06	DT35TB06
to use with AngleBase® & ELLIPTIBase®			

CoCr Abutment			
		30.007	30.054
	Includes screw 19.006 and wax up sleeve for RP and screw19.031 for WP		

All Plastic Castable Abutments			
	Single	28.007	28.054
	Pack of 10	11.007-P10	11.054-P10
		Single unit includes screw 19.006 for RP and screw 19.031 for WP Pack of 10 - screws sold separately	

RPWP

Impression Transfer

		17.007	17.050
		Impression transfer packaged with open and closed tray screw and transfer.	
		<u>for closed tray only</u>	
	5.0mmD	Single	17.507
		Pack of 4	17.507-P4
	8.0mmD	Single	17.807
		Pack of 4	17.807-P4
Hand tightens onto multiunit			

Titanium Temporary Abutments

	Includes screw 19.006 for 24.007 and screw 19.031 for 24.030	24.007	24.030
---	--	--------	--------

Healing abutments

	screws in using nobel unigrip driver	20.100	20.054
---	--------------------------------------	--------	--------

DESSLoc®

	To see the DESSLoc® range please refer to page 53.	88.007	
---	--	--------	--

Scan Bodys

	for desktop scanning * use indicated screw- sold seperately	50.007 * 19.006	50.047 * 19.031
	for intra-oral scanning	52.007	52.047
	on Ti base	53.107	53.154
	on AngleBase	53.107	53.155

ALL DESS MULTIUNITS TAKE THE RP COMPONENTS (4.8mmD)

		IMPLANT PLATFORM	
		NP	RP
	Straight	1,5mm ●	40.007/1 40.008/1
		2,5mm ●	40.007/2 40.008/2
		3,5mm ●	40.007/3 40.008/3
		4,5mm ●	40.008/4
	Angled 17°	2,5mm ◆	41.707/2 41.708/2
		3,5mm ◆	41.707/3 41.708/3
	30°	3,5mm ◆	43.007/3 43.008/3
		4,5mm ◆	43.007/4 43.008/4
with DLC screw			

		IMPLANT PLATFORM				
		NP	RP	WP		
	Straight	1,0mm ●	40.004/1	40.005/1	40.006/1	
		2,0mm ●	40.004/2	40.005/2	40.006/2	
		3,0mm ●	40.004/3	40.005/3	40.006/3	
		4,0mm ●		40.005/4		
		5,0mm ●		40.005/5		
	Angled 17° 	2,0mm ◀	41.704/2	41.705/2		
		3,0mm ◀	41.704/3	41.705/3		
		4,0mm ◀		41.705/4		
		30° 	4,0mm ◀		43.005/4	
			5,0mm ◀		43.005/5	

with DLC screw

ALL DESS MULTIUNITS TAKE THE RP COMPONENTS (4.8mmD) except for only the Branemark WP

IMPLANT PLATFORM

NP RP WP

Nobel Brånemark®



Straight	1,0mm	●	40.001/1	40.002/1	40.003/1
	2,0mm	●	40.001/2	40.002/2	40.003/2
	3,0mm	●	40.001/3	40.002/3	40.003/3
	4,0mm	●		40.002/4	
	5,0mm	●		40.002/5	



Angled 17°	2,0mm	●	41.701/2	41.702/2	
	3,0mm	●	41.701/3	41.702/3	
	4,0mm	●		41.702/4	



30°	4,0mm	●		43.002/4	
	5,0mm	●		43.002/5	

with DLC screw

IMPLANT PLATFORM

3.4 4.1

3i Osseotite External Hex®



Straight	1,0mm	●	40.009/1	40.017/1	
	2,0mm	●	40.009/2	40.017/2	
	3,0mm	●	40.009/3	40.017/3	
	4,0mm	●		40.017/4	
	5,0mm	●		40.017/5	



Angled 17°	2,0mm	●		41.712/2	
	3,0mm	●		41.712/3	
	4,0mm	●		41.712/4	



30°	3,0mm	●		43.012/3	
	4,0mm	●		43.012/4	
	5,0mm	●		43.012/5	

IMPLANT PLATFORM

NC RC

Straumann® Bone Level



Straight	1,5mm	●	40.043/1	40.044/1
	2,5mm	●	40.043/2	40.044/2
	3,5mm	●	40.043/3	40.044/3
	4,5mm	●		40.044/4



Angled 17°	2,5mm	●	41.743/2	41.744/2
	3,5mm	●	41.743/3	41.744/3



30°	3,5mm	●	43.043/3	43.044/3
	4,5mm	●	43.043/4	43.044/4

IMPLANT PLATFORM

RB/WB

Straumann® BLX



Straight	1,5mm	●	40.062/1	
	2,5mm	●	40.062/2	
	3,5mm	●	40.062/3	



Angled 17°	2,5mm	★	41.762/2	
	3,5mm	★	41.762/3	



30°	3,5mm	★	43.062/3	
	4,5mm	★	43.062/4	

ALL DESS MULTIUNITS TAKE THE RP COMPONENTS (4.8mmD)

Neodent Grand Morse®

	Straight		1,5mm	●		40.035/1
			2,5mm	●		40.035/2
			3,5mm	●		40.035/3
	Angled 17°		2,5mm	●		41.035/2
			3,5mm	●		41.035/3
	30°		3,5mm	●		43.035/3
			4,5mm	●		43.035/4

Zimmer Screw-vent®

	Straight		1,0mm	●		40.019/1	40.020/1
			2,0mm	●		40.019/2	40.020/2
			3,0mm	●		40.019/3	40.020/3
			4,0mm	●		40.019/4	
			5,0mm	●		40.019/5	
	Angled 17°		2,5mm	●		41.719/2	41.720/2
			3,5mm	●		41.719/3	41.720/3
	30°		3,5mm	●		43.019/3	43.020/3
			4,5mm	●		43.019/4	43.020/4

Astra EV™

	Straight		1,5mm	●		40.027/1	40.036/1	40.037/1
			2,5mm	●		40.027/2	40.036/2	40.037/2
			3,5mm	●		40.027/3	40.036/3	40.037/3
	Angled 17°		3,5mm	●		41.727/3	41.736/3	41.737/3
			4,5mm	●		41.727/4	41.736/4	41.737/4
	30°		4,5mm	●		43.027/4	43.036/4	43.037/4
			5,5mm	●		43.027/5	43.036/5	43.037/5

Astra Tech Osseospeed™

	Straight		1,5mm	●		40.024/1	40.025/1
			2,5mm	●		40.024/2	40.025/2
			3,5mm	●		40.024/3	40.025/3
	Angled 17°		3,5mm	○		41.724/3	41.725/3
			4,5mm	○		41.724/4	41.725/4
	30°		4,5mm	○		43.024/4	43.025/4
			5,5mm	○		43.024/5	43.025/5

IMPLANT PLATFORM

3.6 4.2 4.8

IMPLANT PLATFORM

AQUA LILAC

3.5-4.0 4.5-5.0

ALL DESS MULTIUNITS TAKE THE RP COMPONENTS (4.8mmD)

		IMPLANT PLATFORM			
		MINI	STANDARD		
	Straight	1,5mm	●	40.072/1	40.073/1
		2,5mm	●	40.072/2	40.073/2
		3,5mm	●	40.072/3	40.073/3
		4,5mm	●		
	Angled 17°	2,5mm	●	41.072/2	41.073/2
		3,5mm	●	41.072/3	41.073/3
	Angled 30°	3,5mm	●	43.072/3	43.073/3
		4,5mm	●	43.072/4	43.073/4

		IMPLANT PLATFORM			
		3.5	4.5		
	Straight	1,0mm	●	40.019/1	40.020/1
		2,0mm	●	40.019/2	40.020/2
		3,0mm	●	40.019/3	40.020/3
		4,0mm	●	40.019/4	
		5,0mm	●	40.019/5	
	Angled 17°	2,5mm	●	41.719/2	41.720/2
		3,5mm	●	41.719/3	41.720/3
	Angled 30°	3,5mm	●	43.019/3	43.720/3
		4,5mm	●	43.019/4	43.720/4

		IMPLANT PLATFORM			
		SP	WP		
	Straight	1,0mm	●	40.019/1	40.020/1
		2,0mm	●	40.019/2	40.020/2
		3,0mm	●	40.019/3	40.020/3
		4,0mm	●	40.019/4	
		5,0mm	●	40.019/5	
	Angled 17°	2,5mm	●	41.719/2	41.720/2
		3,5mm	●	41.719/3	41.720/3
	Angled 30°	3,5mm	●	43.019/3	43.020/3
		4,5mm	●	43.019/4	43.020/4

Dentium Superline™ & Implantium

	Straight	1,2mm	●	40.090/1
		2,5mm	●	40.090/2
		3,5mm	●	40.090/3
	Angled 17°	2,5mm	●	41.790/2
		3,5mm	●	41.790/3
	Angled 30°	2,5mm	●	43.090/2
		3,5mm	●	43.090/3

RP / YELLOW / GREEN / BLUE / RED
3.8 3.6 4.0 4.5 5.0



Titanium Analog - For 3D Printed Models



Single	14.054
Pack of 10	14.054-P10

Our analogs are suitable for both 3D and stone models.

Impression Transfer



17.054

Retention Kits



* divergence up to 40°



divergence between implants			
up to	20°	x1	DLRK-P1
	20°	x2	DLRK-P2
up to	40°	x1	DLER-P1
	40°	x2	DLER-P2

Male Inserts

					divergence between implants	
	up to	20°	x4	1.5 lbs / 680g		DLN2-P4
	up to	20°	x4	3.0 lbs / 1360g		DLN3-P4
	up to	20°	x4	5.0 lbs / 2270g		DLN4-P4
	up to	40°	x4	1.5 lbs / 450g		DLN1-P4
	up to	40°	x4	2.0 lbs / 910g		DLN6-P4
	up to	40°	x4	4.0 lbs / 1810g		DLN7-P4

Tools



DESSLoc Removal & Insertion Tool PHLOC

SPECIAL TOOLS



DESSLoc® DTLOC



Adaptor from Straumann® torque wrench to square connection DTSTSQ



Nobel Biocare™ Adapter from 8mm to square connection DTNBSQ


Nobel Active® / NobelReplace® CC

		NP	RP
		PURPLE 3.5	YELLOW 4.3-5.0
	1,0mm	88.039/1	88.040/1
	2,0mm	88.039/2	88.040/2
	3,0mm	88.039/3	88.040/3
	4,0mm	88.039/4	88.040/4

Nobel Replace Select™

		PURPLE 3.5	YELLOW 4.3
	1,0mm	88.004/1	88.005/1
	2,0mm	88.004/2	88.005/2
	3,0mm	88.004/3	88.005/3
	4,0mm	88.004/4	88.005/4

Nobel Brånemark®

		3.5	4.1
	1,0mm	88.001/1	88.002/1
	2,0mm	88.001/2	88.002/2
	3,0mm	88.001/3	88.002/3
	4,0mm		88.002/4
	5,0mm		88.002/5

3i Osseotite External Hex

		PURPLE 3.4	BLUE 4.1
	1,0mm	88.011/1	88.002/1
	2,0mm	88.011/2	88.002/2
	3,0mm	88.011/3	88.002/3
	4,0mm	88.011/4	88.002/4
	5,0mm		88.002/5

3i Certain Internal Connection

		PURPLE 3.4	BLUE 4.1
	1mm	88.014/1	88.015/1
	2mm	88.014/2	88.015/2
	3mm	88.014/3	88.015/3
	4mm	88.014/4	88.015/4

Straumann® Soft Tissue Level & synOcta®

		4.8 RN
	1,0mm	88.009/1
	2,0mm	88.009/2
	3,0mm	88.009/3
	4,0mm	88.009/4

Straumann® Bone Level

		NC 3.3	RC 4.1
	1,0mm		88.042/1
	2,0mm	88.041/2	88.042/2
	3,0mm	88.041/3	88.042/3
	4,0mm	88.041/4	88.042/4

Zimmer Screw Vent®

		GREEN / 3.5	PURPLE / 4.5
	1,0mm	88.017/1	88.018/1
	2,5mm	88.017/2	88.018/2
	3,5mm	88.017/3	88.018/3
	4,5mm	88.017/4	88.018/4
	5,5mm	88.017/5	

Astra™ EV

		PURPLE 3.6	YELLOW 4.2	BLUE 4.8
	1,0mm	88.059/1	88.060/1	88.061/1
	2,0mm	88.059/2	88.060/2	88.061/2
	3,0mm	88.059/3	88.060/3	88.061/3
	4,0mm	88.059/4	88.060/4	88.061/4

Astra Tech Osseospeed™

		YELOW 3.0	AQUA 3.5-4.0	LILAC 4.5-5.0
	1,0mm	88.023/1	88.024/1	88.025/1
	2,0mm	88.023/2	88.024/2	88.025/2
	3,0mm	88.023/3	88.024/3	88.025/3
	4,0mm	88.023/4	88.024/4	88.025/4
5,0mm				88.025/5

Dentsply friadent® XIVE®

		SILVER 3.4	YELLOW 3.8
	1,0mm	88.036/1	88.037/1
	2,0mm	88.036/2	88.037/2
	3,0mm	88.036/3	88.037/3
	4,0mm	88.036/4	88.037/4

Ankylos® C/X

		
	2,0mm	88.046/2
	3,0mm	88.046/3
	4,0mm	88.046/4

Camlog®

		YELLOW 3.8	PINK 4.3
	1,0mm	88.065/1	88.066/1
	2,0mm	88.065/2	88.066/2
	3,0mm	88.065/3	88.066/3
	4,0mm	88.065/4	88.066/4

BioHorizons® Internal Hex

		YELLOW 3.5	GREEN 4.5
	1,0mm	88.017/1	88.018/1
	2,5mm	88.017/2	88.018/2
	3,5mm	88.017/3	88.018/3
	4,5mm	88.017/4	88.018/4
	5,5mm	88.017/5	

MIS® Seven

		PURPLE / STD	GREEN / WIDE
	1,0mm	88.017/1	88.018/1
	2,5mm	88.017/2	88.018/2
	3,5mm	88.017/3	88.018/3
	4,5mm	88.017/4	88.018/4
	5,5mm	88.017/5	

Multi-Unit Connection®

		MULTI-UNIT 4.8
		88.007

Bar Attachment

		88.100

PRE-MILLED BLANK

All pre milled blank abutments include final screw

		PURPLE YELLOW BLUE				
NobelActive® NobelReplace® CC		3.0	NP	RP	WP	
	Ti	ø10mm	61.068	61.041	61.042	61.069
		ø14mm		62.041	62.042	62.069
	CoCr	ø14mm	71.068	71.041	71.042	71.069

		PURPLE YELLOW BLUE				
NobelReplace Select™		3.5	4.3	5.0	6.0	
	Ti	ø10mm	61.004	61.005	61.006	61.097
		ø14mm		62.005	62.006	62.097
	CoCr	ø14mm	71.004	71.005	71.006	71.097

Nobel Branemark®		3.5	4.1	5.1	
	Ti	ø10mm	61.001	61.002	61.003
	CoCr	ø14mm	71.001	71.002	71.003

			PURPLE	BLUE	YELLOW
3i Osseotite External Hex [®]			3.4	4.1	5.0
	Ti	ø10mm	61.011	61.012	61.013
	CoCr	ø14mm	71.011	71.012	71.013

3i Certain®		PURPLE 3.4	BLUE 4.1	YELLOW 5.0	
	Ti	ø10mm	61.014	61.015	61.016
		ø14mm		62.015	62.016
	CoCr	ø14mm	71.014	71.015	71.016

Straumann® Soft Tissue Level & synOcta®			NNC	RN	WN
			3.5	4.8	6.5
	Ti	ø10mm	61.070	61.009	61.045
		ø14mm		62.009	62.045
	CoCr	ø14mm	71.070	71.009	71.045

Straumann® Bone Level		NC	RC	
		3.3	4.1	
	Ti	ø10mm	61.043	61.044
		ø14mm	62.043	62.044
	CoCr	ø14mm	71.043	71.044

Straumann® BLX		RB	WB
		3.5	4.5
	Ti	ø10mm	61.062
		ø14mm	62.063
	CoCr	ø14mm	71.063

Neodent Grand Morse®

	Ti	ø10mm				61.035
		ø14mm				62.035
	CoCr	ø14mm				71.035

Zimmer Screw-Vent®

			GREEN 3.5	PURPLE 4.5	YELLOW 5.7	
	Ti	ø10mm	61.017	61.018	61.019	
		ø14mm	62.017	62.018	62.019	
	CoCr	ø14mm	71.017	71.018	71.019	

Astra™ EV

			PURPLE 3.6	YELLOW 4.2	BLUE 4.8	
	Ti	ø10mm	61.059	61.060	61.061	
		ø14mm	62.059	62.060	62.061	
	CoCr	ø14mm	71.059	71.060	71.061	

Astra Tech Osseospeed™

			3.0	AQUA 3.5-4.0	LILAC 4.5-5.0	
	Ti	ø10mm	61.023	61.024	61.025	
		ø14mm		62.024	62.025	
	CoCr	ø14mm		71.024	71.025	

Dentsply friadent® XIVE®

			SILVER 3.4	YELLOW 3.8	BLUE 4.5	RED 5.5
	Ti	ø10mm	61.038	61.039	61.040	61.071
		ø14mm		62.039	62.040	62.071
	CoCr	ø14mm	71.038	71.039	71.040	71.071

Ankylos® C/X

	Ti	ø10mm	61.250			
	CoCr	ø14mm				
			71.250			

Hiossen®

				YELLOW MINI	GREEN REGULAR	
	Ti	ø10mm	61.072		61.073	
		ø14mm	62.072		62.073	
	CoCr	ø14mm	71.072		71.073	

Camlog®

			SILVER 3.3	YELLOW 3.8	PINK 4.3	BLUE 5.0
	Ti	ø10mm	61.064	61.065	61.066	61.067
		ø14mm		62.065	62.066	62.067
	CoCr	ø14mm	71.064	71.065	71.066	71.067

Conelog®



		3.3	3.8	5.0
Ti	ø10mm	61.074	61.075	61.076
	ø14mm	62.074	62.075	62.076
	CoCr ø14mm	71.074	71.075	71.076

BioHorizons® Internal Hex



		SILVER 3.0	YELLOW 3.5	GREEN 4.5	BLUE 5.7
Ti	ø10mm	61.099	61.087	61.088	61.089
	ø14mm			62.088	62.089

BioHorizons® External Hex



		YELLOW 3.7	GREEN 4.2	BLUE 5.2
Ti	ø10mm	61.020	61.021	61.022
CoCr	ø14mm		71.021	71.022

MIS® Seven Internal Hex



		PURPLE STD		GREEN WIDE
Ti	ø10mm	61.056	61.095	61.096
CoCr	ø14mm		71.095	71.096

MIS® V3



		BLUE	PURPLE
Ti	ø10mm	61.084	61.082

MIS® C1



		YELLOW	PURPLE	GREEN
Ti	ø10mm	61.056	61.095	61.096

Megagen AnyRidge®



Ti	ø10mm	61.057
CoCr	ø14mm	71.057

Dentium Superline™ & Implantium®



Ti	ø10mm	61.090
	ø14mm	62.090
CoCr	ø14mm	71.090

Blocks




Ti	□15mm L 20mm	63.020
	□15mm L 40mm	63.040
CoCr	□15mm L 20mm	73.020
	□15mm L 40mm	73.040

AngleBase Milling Tool



New



	ØA	ØB	ØC	L1	L2	L3	R
MT3D16L5	3	1,6	0,95	5	51	0,5	0,1
MT4D16L5	4	1,6	0,95	5	51	0,5	0,1
MT6D16L5	6	1,6	0,95	5	51	0,5	0,1

HOLDERS

New

				
HLDSTD	HLDIMES	HLDVHF	HLDVHFR5	HLDYENA
for standard 98mm disc machines	for Imes-Icore® machines	for VHF®	for VHF®	for Yenadent®
Any milli ng machin holder	Coritec 245i touch Coritec 250i touch Coritec 350i / Coritec 350i Loader Coritec 450i Coritec 550i Coritec 650i / Coritec 650i Loader	VHF® S1 VHF® S2 VHF® K5	VHF R5	D6 D15 D43N
				
HLDVHFN4	HLDVHFZ4	HLD3R		
for VHF®	for VHF®	for machines with support/ changer 3R		
VHF® N4	VHF® Z4	D14P DC40 D43N		

CONTACT YOUR VHF PROVIDER FOR DESS HOLDER

Adaptors for ZirkonZahn® holders



HLDZZ

Adaptadores para máquinas ZirkonZahn

M1 abutment
M1 wet heavy metal
M4 wet heavy metal
M5 heavy metal
M6 heavy metal
wet heavy metal



TOOLS

LATCH SCREW DRIVERS

New

HEX	18mm	20mm	24mm	30mm	35mm	compatible with
0,90mm		DTC20H090	DTC24H090			NeodentCM
1,00mm		DTC20H100	DTC24H100		DTC35H100	Ankylos
1,20mm	DTC18H120	DTC20H120	DTC24H120	DTC30H120	DTC35H120	3i Megagen Neodent CM/GM
1,22mm		DTC20H122	DTC24H122		DTC35H122	Xive/Friadent Hiossen
1,25mm		DT20H125	DT4H125			
1,27mm		DTC20H127	DTC24H127	DTC30H127	DTC35H127	Zimmer BioHorizons MIS Astra
1,30mm		DT20S130				Square Driver 1,30mmD

Straumann Drivers- Torx® 06

DT20T06	DT24T06	DT30T06	DT35T06	Straumann
DT20TC06	DT24TC06	DT30TC06	DT35TC06	
DT20TB06	DT24TB06	DT30TB06	DT35TB06	AngleBase

Nobel Unigrip Drivers - UG

DT20U	DT24U	DT30U	DT35U	Nobel
-------	-------	-------	-------	-------

Neodent GM Drivers - NEO

DT24NG	DT30NG	DT35NG	Neodent GM
--------	--------	--------	------------

New

HANDLES AND ADAPTORS

Laboratory handle- for latch drivers	MDTLAA
Thumb Knob handle for Latch Drivers	
7mm	MDTMAN700
9mm	MDTMAN900
Latch style handpiece adaptor for 4x4 Square Connection	
DTANGSQ	
4x4 Square Driver for inserting Straight Multiunit Abutments	
DTMURP	
4x4 Square Driver for latch style driver	
TWADAP	
4x4 Square Driver for DESS LOC	
DTLOC	
Adaptor from Nobel torque wrench to 4x4 square connection	
DTNBSQ	
Adaptor from Straumann torque wrench to 4x4 square connection	
DTSTSQ	
Adaptor from Nobel torque wrench to Latch connection	
DTNBCA	
Adaptor from Straumann torque wrench to Latch connection	
DTSTCA	

Drivers "TC" is a newer version of the "T" driver. For more information on this design please see page 70.



ACTIVE HEX

FITS
NOBELACTIVE® &
NOBELREPLACE® CC



TRI-LOBE

FITS
NOBELREPLACE
SELECT™



CLINICAL SCREWS RECOMMENDED TORQUE

PLATFORM SIZE	DESS REF. N°	DESS RECOMMENDED TORQUE (Ncm)	ORIGINAL EQUIVALENT REF. N°
3.0	19.068	15	Nobel Biocare® 36805
3.0	19.268	15	Nobel Biocare® 37890
3.0	19.668	15	
NP	19.027	35	Nobel Biocare® 31168
NP	19.227	35	Nobel Biocare® 29282
NP	19.101	35	
NP	19.327	35	
NP	19.627	35	
NP AngleBase®	19.441	35	
RP	19.030	35	Nobel Biocare® 28815
RP	19.230	35	Nobel Biocare® 36918
RP	19.328	35	
RP	19.630	35	
RP AngleBase®	19.442	35	
On1®	19.021	35	Nobel Biocare® 38751
On1®	19.221	35	Nobel Biocare® 38750
On1®	19.621	35	
NP on implant for Multi-unit®	19.507	35	
RP on implant for Multi-unit®	19.508	35	

NP	19.004	35	Nobel Biocare® 31170
NP	19.204	35	Nobel Biocare® 36818
NP	19.304	35	
NP	19.604	35	
NP AngleBase®	19.404	35	
RP/WP/6.0	19.005	35	Nobel Biocare® 29293
RP/WP/6.0	19.205	35	Nobel Biocare® 29475
RP/WP/6.0	19.305	35	
RP/WP/6.0	19.605	35	
RP/WP AngleBase®	19.405	35	
NP on implant for Multi-unit®	19.504	35	
RP on implant for Multi-unit®	19.505	35	

EXT HEX UNIVERSAL



EXT HEX USA



INT HEX CLICK



OCTAGON

FITS
STRAUMANN®
SOFT TISSUE
LEVEL & synOcta®



PLATFORM SIZE	DESS REF. N°	DESS RECOMMENDED TORQUE (Ncm)	ORIGINAL EQUIV EQUIVALENT REF. N°
NP	19.027	35	Nobel Biocare® 31168
NP	19.227	35	Nobel Biocare® 29282
NP	19.327	35	
NP	19.101	35	
NP	19.627	35	
NP AngleBase®	19.401	35	
RP	19.028	35	Nobel Biocare® 29290
RP	19.228	35	Nobel Biocare® 29283
RP	19.328	35	
RP	19.628	35	
RP AngleBase®	19.402	35	
WP	19.029	35	Nobel Biocare® 31169
WP	19.229	35	Nobel Biocare® 29284
WP	19.329	35	
WP	19.629	35	
NP on implant for Multi-unit®	19.501	35	
RP on implant for Multi-unit®	19.502	35	

NP/RP/WP	19.034	20	Biomet 3i® UNIHT
NP/RP/WP	19.634	20	
NP/RP/WP AngleBase®	19.402	20	
RP on implant for Multi-unit®	19.512	20	
Low Profile® Abutment ø4.8	19.106	10	Biomet 3i® LPCTSH

NP/RP/WP	19.011	20	Biomet 3i® IUNIHT
NP/RP/WP	19.311	20	
NP/RP/WP	19.611	20	
NP/RP/WP AngleBase®	19.438	20	
Low Profile® Abutment ø4.8	19.106	10	Biomet 3i® LPCTSH

NNC	19.070	35	Straumann® AG 048.314
NNC	19.670	35	
RN/WN on implant	19.026	35	Straumann® AG 048.356
RN/WN on implant	19.102	35	Straumann® AG 049.181
RN/WN on implant	19.328	35	
RN/WN on implant	19.626	35	
RN/WN AngleBase®	19.446	35	
RN/WN on synOcta® abutment	19.009	15	Straumann® AG 048.350
RN/WN on synOcta® abutment	19.103	15	
RN/WN on synOcta® abutment	19.609	15	

CONICAL
BL

FITS
STRAUMANN®
BONE LEVEL



PLATFORM SIZE	DESS REF. N°	DESS RECOMMENDED TORQUE (Ncm)	ORIGINAL EQUIVALENT REF. N°
NC on implant	19.032	35	Straumann® AG 025.2908
NC on implant	19.632	35	
RC on implant	19.033	35	Straumann® AG 025.4908
RC on implant	19.633	35	
NC/RC on implant	19.037	35	Straumann® AG 025.4900
NC/RC on implant	19.104	35	
NC/RC on implant	19.311	20	Straumann® AG 025.4906
NC/RC AngleBase®	19.443	35	
NC/RC on implant for Multi-unit®	19.544	35	Straumann® AG 025.0900-04
NC/RC on Multi-base® abutment	19.041	15	
NC/RC on Multi-base® abutment	19.105	15	Straumann® AG 023.4763
NC/RC on SRA® abutment	19.098	15	

CONICAL BLX

FITS
STRAUMANN®
BLX



RB/WB	19.054	35	
RB/WB	19.700	35	
RB/WB AngleBase®	19.441	35	

NEO GM

FITS
NEODENT® GRAND
MORSE®



GM	19.051	20	
GM AngleBase®	19.444	20	

INT HEX
USA

FITS
ZIMMER
SCREW-VENT®



NP/RP/WP	19.012	30	Zimmer® URS2
NP/RP/WP	19.012L	30	
NP/RP/WP	19.042	30	Zimmer® MHLAS
NP/RP/WP	19.642	30	
NP/RP/WP	19.312	30	Zimmer® MHLAS
NP/RP/WP AngleBase®	19.417	30	
NP (MIS®)	19.327	35	
NP on implant for Multi-unit®	19.519	30	
RP on Multi-unit® abutment ø4.8	19.110	15	
RP on Multi-unit® abutment ø4.8	19.300	15	

CONIC
EVO

FITS
ASTRA TECH
IMPLANT
SYSTEM™ EV



PLATFORM SIZE	DESS REF. N°	DESS RECOMMENDED TORQUE (Ncm)	ORIGINAL EQUIVALENT REF. N°
3.6 on implant	19.059	25	
3.6 on implant	19.659	25	
3.6 AngleBase®	19.441	25	
4.2 on implant	19.060	25	
4.2 on implant	19.660	25	Astra Tech® 25.205
4.2 AngleBase®	19.460	25	
4.8 on implant	19.061	25	
4.8 on implant	19.661	25	
4.8 AngleBase®	19.461	25	
UniAbutment® EV	19.058	15	
UniAbutment® EV	19.658	15	

INT HEX

CONIC

FITS™
ASTRA TECH
OSSEOSPEED



NP on implant	19.015	15	Dentsply® 24189
NP on implant	19.615	15	
RP on implant	19.016	20	Dentsply® 24449
RP on implant	19.616	20	
RP AngleBase®	19.441	20	
WP on implant	19.017	25	Dentsply® 24209
WP on implant	19.617	25	
WP AngleBase®	19.425	25	
UniAbutment® 20°/45°	19.018	15	Dentsply® 22435
UniAbutment® 20°/45°	19.618	15	

INT HEX FD

FITS
DENTSPLY
FRIADENT XIVE®



NP/RP/WP	19.024	24	Dentsply® 454305
NP/RP/WP	19.624	24	
NP/RP/WP AngleBase®	19.438	24	

INT ANK

FITS
DENTSPLY
ANKYLOS® C/X



All platforms	19.048	15	
All platforms	19.050	15	
All platforms	19.648	15	

CONIC
OSS

FITS
OSSTEM® TS



PLATFORM SIZE	DESS REF. N°	DESS RECOMMENDED TORQUE (Ncm)	ORIGINAL EQUIVALENT REF. N°
Mini on implant	19.072	20	Osstem® GSABSST
Mini on implant	19.672	20	
Mini AngleBase®	19.444	20	
Standard on implant	19.073	30	
Standard on implant	19.673	30	
Standard AngleBase®	19.461	30	

INT CAM

FITS
CAMLOG®



3,3 / 3,8 / 4,3 on implant	19.044	20	Camlog® J4005. 1601
3,3 / 3,8 / 4,3 on implant	19.644	20	
3,8 / 4,3 AngleBase®	19.444	20	
5.0 on implant	19.045	20	Camlog® J4005. 2001
5.0 on implant	19.645	20	

INT CON

FITS
CONELOG®



3,3 / 3,8 / 4,3 on implant	19.074	20	Camlog® C4005. 1601
3,3 / 3,8 / 4,3 on implant	19.674	20	
5.0 on implant	19.076	20	Camlog® C4005. 2001
5.0 on implant	19.676	20	

INT HEX
BH

FITS
BIOHORIZONS®
INTERNAL



On implant	19.043	30	BioHorizons® PXAS
On implant	19.043	30	
Multi-unit® RP ø4.8	19.006	15	Nobel Biocare® 29287
Low Profile® abutment ø4.8	19.106	10	Biomet 3i® LPCTSH
Multi-unit® RP ø4.8	19.110	15	
Multi-unit® RP ø4.8	19.206	15	Nobel Biocare® 29285
Multi-unit® RP ø4.8 AngleBase®	19.306	15	
Multi-unit® RP ø4.8	19.606	15	
MIS® Multi-unit® RP ø4.8	19.300	15	
Multi-unit® WP ø6.0	19.031	15	Nobel Biocare® 31163
Multi-unit® WP ø6.0	19.231	15	Nobel Biocare® 29286
Multi-unit® WP ø6.0 AngleBase®	19.331	15	
Multi-unit® WP ø6.0	19.631	15	

EXT HEX
BH

FITS
BIOHORIZONS®
EXTERNAL



PLATFORM SIZE	DESS REF. N°	DESS RECOMMENDED TORQUE (Ncm)	ORIGINAL EQUIVALENT REF. N°
NP	19.013	30	BioHorizons® 130-050
NP	19.312	30	
RP/WP	19.014	30	BioHorizons® 140-050
RP/WP	19.328	35	

INT HEX
MI


FITS
MIS® SEVEN



NP	19.056	30	MIS® MN-S0160
NP	19.656	30	
SP/WP	19.081	30	MIS® MD-S0220
SP/WP	19.681	30	
SP/WP AngleBase®	19.417	30	
Multi-unit® MIS®	19.085	25	
Multi-unit® MIS®	19.685	25	MIS® MU-S0220

CONIC
MC


FITS
MIS® C1



NP	19.056	30	MIS® MN-S0160
NP	19.656	30	
SP/WP	19.081	30	MIS® MD-S0220
SP/WP	19.681	30	
Multi-unit® MIS®	19.085	25	
Multi-unit® MIS®	19.685	25	MIS® MU-S0220

CONIC
MV

FITS
MIS® V3



NP	19.056	30	MIS® MN-S0160
NP	19.656	30	
SP	19.081	30	MIS® MD-S0220
SP	19.681	30	
Multi-unit® MIS®	19.085	25	
Multi-unit® MIS®	19.685	25	MIS® MU-S0220

CONIC
ANYR

FITS
MEGAGEN
AnyRidge®



All platforms	19.057	35	Megagen® AANMSF
All platforms	19.657	35	
All platforms AngleBase®	19.460	35	

CONIC DENT





On implant	19.090	30	Dentium® ASC 2045
On implant	19.690	30	
On abutment	19.066	10	Dentium® SRS 18T
On abutment	19.666	10	

MULTI-UNIT





NP	19.006	15
NP	19.206	15
NP	19.606	15
WP	19.031	15
WP	19.231	15
WP	19.631	15
RP AngleBase®	19.306	15
WP AngleBase®	19.631	15

DESS USA SCREW SELECTION GUIDE



BRAND	IMPLANT SYSTEM	DESS CONNECTION	REF	DESCRIPTION	HEAD DIAMETER	HEAD BASE ANGLE	OVERALL LENGTH	HEAD LENGTH	DISTANCE FROM HEAD TO THREAD	THREAD LENGTH	THREAD KIND	TORQUE	DRIVER SOCKET
Nobel Biocare®	Brånemark®	External Hex Universal	19.027	Screw UG NP	2,23	45°	7,35	2,00	2,60	2,75	M1,6	35 N/cm	Unigrip®
			19.028	Screw UG RP	2,53	0°	7,35	1,42	3,12	2,81	M2,0	35 N/cm	
			19.029	Screw UG WP	2,85	0°	7,10	1,40	2,70	3,00	M2,5	35 N/cm	
			19.101	Screw UG NP flat seat	2,23	0°	7,35	1,90	2,70	2,75	M1,6	35 N/cm	
			19.227	DLC Screw UG NP	2,23	45°	7,35	2,00	2,60	2,75	M1,6	35 N/cm	
			19.228	DLC Screw UG RP	2,53	0°	7,35	1,42	3,12	2,81	M2,0	35 N/cm	
			19.229	DLC Screw UG WP	2,85	0°	7,10	1,40	2,70	3,00	M2,5	35 N/cm	
	Replace® Select	Tri-lobed	19.004	Screw UG NP	2,45	0°	8,30	1,70	4,00	2,60	M1,8	35 N/cm	
			19.005	Screw UG RP and WP	2,45	0°	9,90	1,50	4,90	3,50	M2,0	35 N/cm	
			19.204	DLC Screw UG NP	2,45	0°	8,30	1,70	4,00	2,60	M1,8	35 N/cm	
	Active® / Replace® Conical Connection	Active Hex	19.205	DLC Screw UG RP and WP	2,45	0°	9,90	1,50	4,90	3,50	M2,0	35 N/cm	
			19.027	Screw UG NP	2,23	45°	7,35	2,00	2,60	2,75	M1,6	35 N/cm	
			19.030	Screw UG RP	2,53	54°	7,35	1,57	2,97	2,81	M2,0	35 N/cm	
			19.068	Screw UG 3.0	1,82	55°	6,72	1,82	2,45	2,45	M1,4	15 N/cm	
			19.227	DLC Screw UG NP	2,23	45°	7,35	2,00	2,60	2,75	M1,6	35 N/cm	
			19.230	DLC Screw UG RP	2,53	54°	7,35	1,57	2,97	2,81	M2,0	35 N/cm	
Straumann®	Tissue Level	Octagon	19.009	Screw "Torx®" on abutment	2,50	45°	4,50	1,20	0,00	3,30	M2,0	15 N/cm	TORX® 06
			19.026	Screw "Torx®"	2,58	75°	6,70	3,20	1,68	2,00	M2,0	35 N/cm	
			19.070	Screw "Torx®" NNC	2,10	75°	7,10	1,85	3,35	1,70	M1,6	35 N/cm	
			19.102	Screw "Torx®" flat seat	2,40	0°	7,50	1,55	3,25	2,70	M2,0	35 N/cm	
			19.103	Screw "Torx®" flat seat (on abutment)	2,50	0°	5,10	1,55	0,00	3,55	M2,0	15 N/cm	
	Bone Level	Conical BL	19.032	Screw Torx NP	2,20	75°	7,85	2,31	3,84	1,70	M1,6	35 N/cm	
			19.033	Screw Torx RP	2,20	75°	7,85	2,31	3,84	1,70	M1,6	35 N/cm	
			19.037	Screw Torx NC/RC	2,20	75°	7,85	2,31	3,84	1,70	M1,6	35 N/cm	
			19.041	Screw "Torx®" on abutment	2,20	45°	4,50	1,60	0,30	2,60	M1,4	15 N/cm	
			19.098	Screw for SRA	2,20	55°	3,65	1,65	0,20	1,80	M1,4	15 N/cm	
			19.104	Screw "Torx®" flat seat	2,20	0°	8,15	1,55	4,90	1,70	M1,6	35 N/cm	
			19.105	Screw "Torx®" flat seat (on abutment)	2,20	0°	4,50	1,50	0,00	3,00	M1,4	15 N/cm	
			19.015	Screw hexagonal 1,27mm NP (3,0)	2,10	30°	7,50	2,00	3,00	2,50	M1,4	15 N/cm	
			19.016	Screw hexagonal 1,27mm RP	2,35	0°	8,25	2,00	3,22	3,02	M1,6	20 N/cm	
Dentsply Sirona®	Ossespeed™	Internal Hex Conic	19.017	Screw hexagonal 1,27mm WP	2,35	0°	8,30	2,00	4,80	3,50	M2,0	25 N/cm	HEX .050 1,27
			19.018	Screw on abutment hexagonal 1,27mm	2,30	0°	4,80	1,90	0,00	2,90	M1,4	15 N/cm	
			19.058	Screw hex. 1,27mm on Uniabutment EV	2,35	18°	4,49	1,71	0,00	2,50	M1,8	15 N/cm	
	Astra EV	Conic EVO	19.059	Screw hexagonal 1,27mm EV 3.6	2,15	60°	7,46	2,12	1,34	4,00	M1,6	25 N/cm	
			19.060	Screw hexagonal 1,27mm EV 4.2	2,35	40°	7,64	1,85	1,79	4,00	M1,8	25 N/cm	
			19.061	Screw hexagonal 1,27mm EV 4.8	2,35	0°	7,60	1,60	2,00	4,00	M2,0	25 N/cm	
	Xive®	Internal Hex FD	19.024	Screw hexagonal 1,22mm	2,10	45°	9,80	3,06	3,29	3,45	M1,6	25 N/cm	HEX .048 1,22mm
			19.048	Screw hexagonal 1,00mm	2,35	60°	7,87	2,37	2,40	3,10	M1,8	15 N/cm	
	Ankylos® C/X	Internal ANK	19.050	Screw hexagonal 1,00mm	2,35	60°	9,37	2,02	3,95	3,40	M1,8	15 N/cm	Hex 1,00mm
Zimmer -Biomet®	3i™ Osseotite®	External Hex USA	19.034	Screw hexagonal 1,20mm	2,55	0°	7,85	1,65	0,70	5,50	M2,0	20 N/cm	HEX .048 1,20mm
	3i™ Certain®	Internal Hex "Click"	19.011	Screw hexagonal 1,20mm	2,00	0°	8,40	1,90	3,80	2,70	M1,6	20 N/cm	
	Screw-Vent®	Internal Hex USA	19.012	Screw hexagonal 1,27mm	2,30	19°	8,20	2,55	2,98	2,60	UNF 1-72	30 N/cm	HEX .050 1,27
			19.042	Screw hexagonal 1,27mm	1,95	19°	8,20	2,49	2,95	2,76	UNF 1-72	30 N/cm	

DESS USA SCREW SELECTION GUIDE



BRAND	IMPLANT SYSTEM	DESS CONNECTION	REF	DESCRIPTION	HEAD DIAMETER	HEAD BASE ANGLE	OVERALL LENGTH	HEAD LENGTH	DISTANCE FROM HEAD TO THREAD	THREAD LENGTH	THREAD KIND	TORQUE	DRIVER SOCKET
MIS®	Seven®	Internal Hex MI	19.056	Screw hexagonal 1,27mm 3.0	2,10	30°	7,63	2,40	2,53	2,70	M1,6	30 N/cm	HEX .050 1.27
			19.081	Screw hexagonal 1,27mm SP/WP	2,20	28°	7,48	2,41	2,82	2,25	UNF 1-72	30 N/cm	
	C1 / V3	Conic MC / MV	19.056	Screw hexagonal 1,27mm NP	2,10	30°	7,63	2,40	2,53	2,70	M1,6	30 N/cm	
			19.081	Screw hexagonal 1,27mm SP/WP	2,20	28°	7,48	2,41	2,82	2,25	UNF 1-72	30 N/cm	
			19.085	Screw hexagonal 1,27mm Multiunit 4.8	2,15	0°	3,50	1,30	0,00	2,20	M1,6	25 N/cm	
BioHorizons®	External	External Hex BH	19.656	Screw hexagonal 1,27mm NP with TIN	2,10	30°	7,63	2,40	2,53	2,70	M1,6	30 N/cm	HEX .050 1.27
			19.013	Screw hexagonal 1,27mm NP	2,50	0°	7,00	1,00	2,10	3,90	UNF 1-72	30 N/cm	
	Internal	Internal Hex BH	19.014	Screw hexagonal 1,27mm RP/WP	2,65	0°	7,10	1,00	2,10	4,00	M2,0	30 N/cm	
			19.043	Screw hexagonal 1,27mm	2,40	45°	7,24	2,42	2,36	2,45	UNF 1-72	30 N/cm	
Camlog®	Camlog® system	Internal CAM	19.044	Screw hexagonal 1,27mm NP	2,10	45°	10,45	2,35	4,60	3,50	M1,6	20 N/cm	HEX .050 1.27
			19.045	Screw hexagonal 1,27mm WP	2,50	45°	10,45	2,35	4,60	3,50	M2,0	20 N/cm	
	Conelog® system	Internal CON	19.074	Screw hexagonal 1,27mm 3.3, 3.8, 4.3	2,26	30°	8,87	2,08	3,99	2,80	M1,6	20 N/cm	
Megagen®	AnyRidge®	Conic ANYR	19.076	Screw hexagonal 1,27mm 5.0	2,67	30°	8,87	2,08	3,99	2,80	M2,0	20 N/cm	HEX .048
			19.057	Screw hexagonal 1,20mm	2,10	0°	8,00	2,00	2,78	3,22	M1,8	35 N/cm	
Dentium®	Superline Implantium®	Conic Dent	19.066	Screw hexagonal 1,28mm MUA	2,30	45°	5,00	1,93	0,57	2,50	M1,8	10 N/cm	HEX .050
			19.090	Screw hexagonal 1,28mm S-Line	2,30	33°	10,20	2,05	4,85	3,30	M2,0	30 N/cm	
Hiossen Implant®	ETI/ TS®	Conic OSS	19.072	Screw hexagonal 1,22mm Mini	2,22	45°	10,18	2,67	5,41	2,10	M1,6	20 N/cm	HEX .048
			19.073	Screw hexagonal 1,22mm Standard	2,33	45°	8,35	2,60	2,95	2,80	M2,0	30 N/cm	
Multiunit		Multiunit	19.006	Screw UG on Multi-Unit® RP	2,00	0°	3,45	1,35	0,00	2,10	M1,4	15 N/cm	Unigrip®
			19.031	Screw UG on Multi-Unit® WP	2,50	0°	4,35	1,65	0,00	2,70	M1,8	15 N/cm	Unigrip®
			19.106	Screw H1.20 on 3i® Munit® RP	2,25	0°	4,48	1,75	0,53	2,20	M1,4	10 N/cm	HEX .048
			19.110	Screw H1.27mm on Multi-Unit® RP	2,00	0°	3,45	1,35	0,00	2,10	M1,4	15 N/cm	HEX .050
			19.206	DLC Screw UG on Multi-Unit® RP	2,00	0°	3,45	1,35	0,00	2,10	M1,4	15 N/cm	Unigrip®
		Screws for Angled MUA	19.231	DLC Screw UG on Multi-Unit® WP	2,50	0°	4,35	1,65	0,00	2,70	M1,8	15 N/cm	Unigrip®
			19.501	DLC screw replacement for ref. 41.701/x	2,45	0°	6,00	0,60	2,70	2,70	M1,6	20 N/cm	Unigrip®
			19.502	DLC screw for ref. 41.702/x & 43.002/x	2,45	0°	6,00	0,60	2,90	2,50	M2,0	20 N/cm	Unigrip®
			19.504	DLC screw replacement for ref. 41.704/x	2,45	0°	7,15	0,60	3,90	2,65	M1,8	35 N/cm	Unigrip®
			19.505	DLC screw for ref. 41.705/x & 43.005/x	2,45	0°	9,00	0,60	4,90	3,50	M2,0	20 N/cm	Unigrip®
			19.507	DLC screw for ref. 41.707/x & 43.007/x	2,45	0°	6,70	0,60	2,78	3,32	M1,6	35 N/cm	Unigrip®
			19.508	DLC screw for ref. 41.708/x & 43.008/x	2,45	0°	6,40	0,60	3,30	2,50	M2,0	20 N/cm	Unigrip®
			19.519	Screw for ref. 41.719/x & 43.019/x	2,45	0°	6,00	0,60	2,80	2,60	UNF 1-72	30 N/cm	HEX .050
			19.544	Screw ref. 41.743-44/x & 43.043-44/x	2,45	0°	7,20	0,60	2,85	3,75	M1,6	35 N/cm	TORX® 06
			19.401	Screw for DESS AngleBase M1,6	2,10	0°	6,10	1,40	2,50	2,20	M1,6	35 N/cm	TORX® 06 Ball Angle Screw Channel Driver
			19.402	Screw for DESS AngleBase M2,0	2,50	0°	5,90	1,40	1,69	2,81	M2,0	35 N/cm	
			19.404	Screw for DESS AngleBase M1,8	2,30	0°	7,20	1,40	3,55	2,25	M1,8	35 N/cm	
			19.405	Screw for DESS AngleBase M2,0	2,45	0°	8,05	1,40	4,60	2,05	M2,0	35 N/cm	
			19.417	Screw for DESS AngleBase UNF-1-72	2,30	0°	4,90	1,40	2,90	2,00	UNF 1-72	30 N/cm	
			19.425	Screw for DESS AngleBase M2,0	2,45	0°	8,90	1,40	5,45	2,05	M2,0	25 N/cm	
			19.438	Screw for DESS AngleBase M1,6	2,10	0°	7,90	1,40	4,50	2,00	M1,6	25 N/cm	
			19.441	Screw for DESS AngleBase M1,6	2,10	0°	7,30	1,40	3,70	2,20	M1,6	35 N/cm	
			19.442	Screw for DESS AngleBase M2,0	2,45	0°	7,20	1,40	3,75	2,05	M2,0	35 N/cm	
			19.443	Screw for DESS AngleBase M1,6	2,10	0°	8,15	1,40	4,75	2,00	M1,6	35 N/cm	
			19.444	Screw for DESS AngleBase M1,6	2,10	0°	9,20	1,40	5,60	2,20	M1,6	20 N/cm	
			19.446	Screw for DESS AngleBase M2,0	2,45	0°	6,70	1,40	3,25	2,05	M2,0	35 N/cm	
			19.460	Screw for DESS AngleBase M1,8	2,30	0°	7,42	1,40	4,02	2,00	M1,8	25 N/cm	
			19.461	Screw for DESS AngleBase M2,0	2,45	0°	7,75	1,40	4,30	2,05	M2,0	25 N/cm	



- non-engaging
- engaging

CEREC COMPATIBLE TIBASE SOLUTION

QUICK REFERENCE GUIDE

			PRODUCT REFERENCES								
CONNECTION	SOFTWARE	GH1		GH2		GH3		SCREW	TORQUE	PLATFORM	
											
NOBELACTIVE® & NOBELREPLACE® CC	L	NB A 4.5 L	58.041	59.041			58.041/3	59.041/3	19.027	35 Ncm	NP
		NB A 5.0 L	58.042	59.042			58.042/3	59.042/3	19.030		RP
NOBELREPLACE SELECT™	L	NB RS 3.5 L	58.004	59.004	58.004/2	59.004/2	58.004/3	59.004/3	19.004	35 Ncm	NP (3.5)
		NB RS 4.3 L	58.005	59.005	58.005/2	59.005/2	58.005/3	59.005/3	19.005		RP (4.3)
		NB RS 5.0 L	58.006	59.006	58.006/2	59.006/2	58.006/3	59.006/3	19.005		WP (5.0)
		NB RS 6.0 L	58.097	59.097	58.097/2	59.097/2	58.097/3	59.097/3	19.005		6.0
NOBEL BRÄNEMARK®	L	NB B 3.4 L		59.001	58.001/2	59.001/2	58.001/3	59.001/3	19.027	35 Ncm	NP (3.4)
		NB B 4.1 L		59.002	58.002/2	59.002/2	58.002/3	59.002/3	19.028		RP (4.1)
3i OSSEOTITE™	L	B 0 3.4 L		59.011	58.011/2	59.011/2	58.011/3	59.011/3	19.034	20 Ncm	3.4
		B 0 4.1 L		59.012	58.012/2	59.012/2	58.012/3	59.012/3	19.034		4.1
		B 0 5.0 L		59.013	58.013/2	59.013/2	58.013/3	59.013/3	19.034		5.0
3i CERTAIN®	S	B C 3.4 S	58.014	59.014	58.014/2	59.014/2	58.014/3	59.014/3	19.011	20 Ncm	3.4
	L	B C 4.1 L	58.015	59.015	58.015/2	59.015/2	58.015/3	59.015/3	19.011		4.1
		B C 5.0 L	58.016	59.016	58.016/2	59.016/2	58.016/3	59.016/3	19.011		5.0
STRAUMANN® SOFT TISSUE LEVEL & synOCTA®		S SO 4.8 L	58.009	59.009					19.026		RN
		S SO 6.5 L	58.045	59.045					19.026		WN
STRAUMANN® BONE LEVEL	L	S BL 3.3 L	58.043	59.043			58.043/3	59.043/3	19.032	35 Ncm	NC
		S BL 4.1 L	58.044	59.044			58.044/3	59.044/3	19.033		RC

QUICK REFERENCE GUIDE

			PRODUCT REFERENCES								
			GH1		GH2		GH3				
CONNECTION	SOFTWARE								SCREW	TORQUE	PLATFORM
ZIMMER® SCREW-VENT®	L	Z TSV 3.5 L	58.017	59.017	58.017/2	59.017/2	58.017/3	59.017/3	19.042	30 Ncm	3.5
		Z TSV 4.5 L	58.018	59.018	58.018/2	59.018/2	58.018/3	59.018/3	19.042		4.5
		Z TSV 5.7 L	58.019	59.019	58.019/2	59.019/2	58.019/3	59.019/3	19.042		5.7
ASTRA TECH IMPLANT SYSTEM™ EV	S	AT EV 3.6 GH1 S	58.059	59.059			58.059/3	59.059/3	19.059	25 Ncm	3.6
	L	AT EV 4.2 GH1 L	58.060	59.060			58.060/3	59.060/3	19.060		4.2
		AT EV 4.8 GH1 L	58.061	59.061			58.061/3	59.061/3	19.061		4.8
ASTRA TECH OSSEOSPEED™	L	AT OS 3.5/4.0 L	58.024	59.024			58.024/3	59.024/3	19.016	20 Ncm	3.5/4.0
		AT OS 4.5/5.0 L	58.025	59.025			58.025/3	59.025/3	19.017	25 Ncm	4.5/5.0
DENTSPLY FRIADENT XIVE®	S	FX 3.4 S	58.038	59.038	58.038/2	59.038/2	58.038/3	59.038/3	19.024	24 Ncm	3.4
		FX 3.8 S	58.039	59.039	58.039/2	59.039/2	58.039/3	59.039/3	19.024		3.8
	L	FX 4.5 L	58.040	59.040	58.040/2	59.040/2	58.040/3	59.040/3	19.024		4.5
		FX 5.5 L	58.071	59.071	58.071/2	59.071/2	58.071/3	59.071/3	19.024		5.5
DENTSPLY ANKYLOS™ C/X	S	ANKN C/ GH1 S	58.050	59.050					19.046	15 Ncm	C/X
		ANKN C/ GH2 S			58.050/2	59.050/2			19.046		
OSSTEM® TS	L	O TS 3.5 L	58.072	59.072			58.072/3	59.072/3	19.072	20 Ncm	MINI
		O TS 4.0 L	58.073	59.073			58.073/3	59.073/3	19.0737	30 Ncm	STANDARD
CONELOG®	S	CONELOG 3.3/x	58.074	59.074	58.074/2	59.074/2			19.052	20 Ncm	3.3
		CONELOG 3.8/x	58.075	59.075	58.075/2	59.075/2			19.052		3.8/4.3
	L	CONELOG 5.0/x	58.076	59.076	58.076/2	59.076/2			19.053		5.0
BIOHORIZONS™ INTERNAL	L	BH 3.5 L	58.087	59.087			58.087/3	59.087/3	19.043	30 Ncm	3.5
		BH 4.5 L	58.088	59.088			58.088/3	59.088/3	19.043		4.5
		BH 5.7 L	58.089	59.089			58.089/3	59.089/3	19.043		5.7

Analogs -
Since 2013 all
DESS analogs can
be used for both
stone and digital
models. One is all
you need. Coming
soon is a digital
analog with screw
at the bottom.

68



Our analogs have been carefully designed to match the implant connection to perfection. Manufactured either in stainless steel AISI 303 and Titanium Grade 5 ELI.

- Color coded to match implant platform
- Precise connection guaranteed
- For 3d and traditional clay models

Coming soon

A 3D analog with screw at the bottom to guarantee that the analog reaches the lower end of the model.




DESSLoc®

Superior resistance, better results

Coated with a Zirconium nitride (ZrN), the DESS® overdenture anchor system offers you better resistance with less plaque adhesion improving results and reducing inflammation

- 100% compatible with the most popular system
- Reduces plaque adhesion
- Decreases inflammation

DESSLoc® coated with Zirconium nitride:



aMMP-8 values of fluid from peri-implant sulcus in titanium (Ti) and Zirconium Nitride (ZrN) coated attachments at different times after prosthetic restoration. (Matrix Metalloproteinase-8 Levels in Peri-implant Sulcus Fluid Adjacent to Titanium and Zirconium Nitride Surface, Richard Meissen, MSc1 Mariana Mintcheva, MSc1 Lutz Netuschil, PhD Dr2, Int J Periodontics Restorative Dent 2014; 34:91-95. doi: 10.11607/prd.1504)

Zrn coating is a patent from Clinical House Dental GmbH




SUPERIOR LATCH DRIVERS

Patented design


Parallel Screwdrivers:

- Clamping
- + Durability
- + Screw head protection



Conical Screwdrivers:


- + Clamping
- Durability
- Screw head protection



By taking the best of both designs we created an overall better screwdriver.

DESS® Screwdrivers:

- + Better clamping
- + Durability
- + Screw head protection



All DESS screwdrivers are manufactured in Stainless Hardened Steel for increased durability and wear resistance.



Ti-Bases

Manufactured in Titanium Grade V ELI, and available in different gingival heights, the DESS Ti-Bases are a great alternative to the implant branded products. The main advantages are:



- **Vertical Walls for Maximum Bond Retention**

Regardless of the chimney height a 5 degree taper has a 68% greater chance of debonding because of the tapered wall...all DESS tibases are vertical walls .

- **SelectGrip® enhanced retention**


Patented surface treatment that increases retention.

- **Different Gingival Heights**

Available in different heights to cater for every case.

- **Anatomic profile**

DESS® Ti-bases geometry design has an anatomic profile for better, deeper and seamless restorations.




SelectGrip®

Our patented surface enhancement treatment improves retention up to 500%.

SelectGrip® is a special sandblasted treatment which provides a considerable improvement in retention levels of cemented prosthesis. Applied for both abutments and Ti-bases, the adhesion to the final cemented prosthesis improves significantly.

This graph shows the average pressure supported before the joint between the ZrN sleeve and titanium tibase breaks




*Comparison between two different sets of titanium abutments and titanium bases cemented to a standard ZrN, sleeve. The cement used for the test was Multilink Hybrid Abutment from Ivoclar Vivadent®.

The importance of a factory treatment:

In house sandblasting, not only means more chair time but if its not done properly (bellow the collar) this can increase the chances of gingivitis.



SelectGrip®



In house

AngleBase®

Angular freedom


AngleBase® is our answer to clinical situations where the long-axis of the implant results in an unfavorable location of the prosthetic screw access hole.

- Full recommended torque at maximum angle
- Up to 25° angle.
- 360° freedom positioning.
- Flat seat screw design.
- 621N Bonding retention.

360°

Get the
angle
you
need

The specially designed Torx® based screw and driver concept will retain the ability to use full recommended torque even at full angulation.




ELLIPTIBase


The Revolutionary Ti-Base

Due to the limited interdental space in the anterior teeth section, restoration in these positions has always been a challenge.

Traditional Ti-Bases have to use less material and in consequence a weaker restoration is achieved. ELLIPTIBase® is our solution specially designed for narrow implant platforms.



The 22% reduction in size, achieved on the design of the ELLIPTIBase® allows more space for the restoration whilst keeping the same cementation area.



SelectGrip® patented enhancement surface treatment for up to 500% more retention



RE-ENGINEERED TO PERFECTION

No need for introductions, the C-Base® is our take on the Sirona® Cerec® Ti-base. Based on the pure switch concept we took on the already flawless Sirona Ti-base and enhanced some of its properties.



Different gingival heights

To cater for different clinical situations



SelectGrip® surface treatment

Better bonding properties, reducing up to 15 min of preparation time per unit

Gold anodized coating

Improved finish for better overall aesthetics

Engaging or non engaging

Non engaging, ideal for bridges



Final screw included

No added costs; Pure Switch same as implant manufactured screw

YOUR DESIGN OUR PRECISION


Mill your design with the precision of the DESS® connections. Our range is available for the most popular implant connections in the world.



THE DEFINITIVE MULTI-UNIT



Straight 17° 30°




Multiple
gingival heights



Full shaped connection
provides more stability.

The DESS® Ti Base for
Multi-Unit® is available in two
different heights: 6mm and 8mm.

Full DIY control



The DESS® AngleBase Milling tool will allow you to create screw retained bars or bridges and multiple with an angulated channel using the AngleBase® screw. These restorations are part of the DESS® “Do it Yourself” concept that allow you to design your pieces with the already known precision of the DESS® connection.

The AngleBase Milling Tool replicates the AngleBase® screw seat.

Available in Titanium and Chrome Cobalt in 15mm and lengths of 20mm or 40mm

 ACTIVE HEX <i>FITS</i> NOBELACTIVE® & NOBELREPLACE® CC	 OCTAGON <i>FITS</i> STRAUMANN® SOFT TISSUE LEVEL & synOcta®	 CONIC EVO <i>FITS</i> ASTRA TECH IMPLANT SYSTEM™ EV	 INTERNAL CON <i>FITS</i> CONELOG®	 CONIC IC <i>FITS</i> MEDENTIS ICX®
 TRI-LOBE <i>FITS</i> NOBELREPLACE SELECT™	 CONICAL BL <i>FITS</i> STRAUMANN® BONE LEVEL	 INTERNAL HEX CONIC <i>FITS</i> ASTRA TECH OSSEOSPEED™	 INTERNAL HEX BH <i>FITS</i> BIOHORIZONS® INTERNAL	
 EXTERNAL HEX UNIVERSAL <i>FITS</i> NOBEL BRANEMARK®	 CONICAL BLX <i>FITS</i> STRAUMANN® BLX	 INTERNAL HEX FD <i>FITS</i> DENTSPLY FRIADENT XIVE®	 EXTERNAL HEX BH <i>FITS</i> BIOHORIZONS® EXTERNAL	
 EXTERNAL HEX USA <i>FITS</i> 3i OSSEOTITE®	 NEO GM <i>FITS</i> NEODENT® GRAND MORSE	 CONIC OSS <i>FITS</i> OSSTEM® TS	 INTERNAL HEX MI <i>FITS</i> MIS® SEVEN	
 INTERNAL HEX “CLICK” <i>FITS</i> 3i CERTAIN	 INTERNAL HEX USA <i>FITS</i> ZIMMER SCREW-VENT®	 INTERNAL CAM <i>FITS</i> CAMLOG®	 CONIC ANYR <i>FITS</i> MEGAGEN AnyRidge®	

LIFETIME WARRANTY

Here at DESS we believe in excellence, we put our passion into every single detail of our products.

Since early days we understood that saying this wasn't enough so we pioneered introducing the lifetime warranty in all our products.

We're so confident about the quality and performance of our products that we give as standard a Lifetime warranty* on all the products we manufacture.

That's the DESS commitment.

If you want to know more about our warranty and the 3rd party implants warranty visit our website www.dess-usa.com and download the warranty.

* Subject to conditions, please visit www.dess-usa.com to download our warranty




The Pure Switch which concept is the philosophy upon our products are conceived.

The best example of this are our Screws which have been manufactured using the same dimensions and characteristics as the implant branded screw.

As the key component of every ensemble, DESS® screws have a perfect fit and can be interchangeable with the branded version without any compromise.

Reasons to switch to DESS®

- **One stop shop**
31 compatibilities with the most popular implant
connections
- **European Quality**
ISO 9001, ISO 13485, CE, FDA, Health Canada
- **Lifetime Warranty**
All DESS® manufactured products have a
lifetime warranty*



“Staying true to our history and values empowers our brand”

DESS Abutments company history:

- 1947 Founded in Barcelona. Manufacture of precision mechanical components.
- 1970 70% of production using Swiss-style automatic lathes.
- 1988 First purchase of a Japanese numerical control Okuma CNC lathe.
- 1997 First Japanese Okuma CNC machining.
- 1999 First incorporation of Swiss DECO lathes specialized in micro precision machining.
- 2000 ISO 9001 Certification.
- 2003 100% of production using numerical control CNC machines.
- 2004 First production of dental implant components.
- 2008 Creation of the sanitary department with its own clean room.
 - Sanitary certificate ISO 13485 and CE.
 - Creation of the DESS Dental Smart Solutions brand.
- 2009 SelectGrip® Patent for surface finishes. 500% proven increased bonding retention.
- 2011 DESS already has a presence in 16 countries around the world.
- 2012 Creation of DESS GmbH in Düsseldorf, Germany.
- 2013 DESS already has a presence in 26 countries around the world.
- 2014 DESS LAB starts offering one piece custom abutments.
- 2015 Health Canada approval.
- 2016 We completed the full assortment of angulated Multi-units.
 - We launched a revolutionary intraoral scan abutment, without hole on top.
 - STA® is accepted by the regulatory bodies as the unique surface offering real soft tissue attachment and finally available in Europe.
 - DESS is present in 34 countries worldwide.
- 2017 Achievement of FDA.
 - Relocation of manufacturing plant and offices with a total of 6500m².
- 2018 We launch our angulated screw channel solution “AngleBase®”.
 - Resellers in 41 countries.
 - 100% Increase of the production.
- 2019 Creation of the DiY range.
 - Subsidiary is created in US.
 - Two new compatibilities launched.
- 2020 The ElliptiBase® is born, revolutionary solution for narrow interdental spaces. Three new compatibilities launched.
 - C-Base®, Neodent, Straumann BLX

Since our beginnings quality has been quintessential to our philosophy and is reflected in every product that we manufacture. Our choice of the best Swiss and Japanese milling machines allow us to create the most precise parts. Our quality standards entail measuring to a few microns of precision.



